

MEMORANDUM

TO: Regional Medical Laboratory Clients

FROM: William Fitter, MD, Chief of Clinical Chemistry

Nikki Garner, MS HCA, Manager Processing/Specimen Resolution/Reference Lab

Shelli Lamb, MLT(ASCP), Reference Lab Specialist

DATE: April 21, 2015

SUBJECT: Chromogranin A

Change in methodology, reference range and availability of crossover studies

Effective May 26, 2015, there will be a change in methodology and reporting of the Chromogranin A assay. Due to these changes, re-baseline testing will be offered at no charge from May 26, 2015 to July 24, 2015. The new method will provide a broader reference range and faster turn-around-time.

RE-BASELINE Ordering Information:

For re-baseline studies, please order the following test only: **CHROMGR RB** (**3803559**) There will be no need to order separate tests. Both results will be reported through your normal charting methods.

The specimen requirement for re-baseline studies is as follows:

Two (2) - 1mL serum specimens. Separate serum from cells and refrigerate as soon as possible.

A report can be generated to list all your patients' re-baseline results by contacting RML Client Services at 918-744-2500 or your Sales Representative.

After July 27, 2015, please order the following test:

New Order Mnemonic: CHROMOG A (3803557).

New methodology: Quantitative Enzyme Immunoassay

Reference Range: 0-95 ng/mL

Specimen Required: 1mL (0.5mL min) serum from either an SST or Red no-gel tube.

Separate serum from cells and refrigerate as soon as possible.

Specimen Stability: Room temp=48 hours; Refrigerated=2 weeks; Frozen=6 weeks

CPT code: 86316

If you have questions or concerns, please feel free to contact either William Fitter, MD, Chief of Clinical Chemistry and Toxicology; Nikki Garner, Manager Processing/Specimen Resolution/Reference Lab; or Shelli Lamb, MLT(ASCP), Reference Lab Specialist at 918-744-2500.