

NEW TESTS - Please update your EMR catalog with those appropriate to your practice

New Test Effective : Immediately									
Test Code	Test Name	Mnemonic	Category/Type						
6907557	SARS-CoV-2 by PCR	SARS-CoV2	Group	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
CPT: U0002	Collect: Nasopharyngeal swab AND Oropharyngeal swab. Place both swabs in one collection tube. Specimen Prep: Place in viral transport medium (VTM). Place each specimen in an individually sealed bag. Unacceptable: Specimens not in VTM, Ambient specimens. Stability: Ambient - Unacceptable; Refrigerated - 4 days; Frozen - 1 month Method: Qualitative PCR Set-up: Daily TAT: 1-2 days after test set-up		Prompt	N/A	6907559	SARS-CoV2 PCR	SARS-CoV-2 by PCR	By report	Alpha
				31208-2	6907561	CoV2 Source	SARS-CoV-2 Source	By report	By Report

New Test Effective : Immediately										
Test Code	Test Name	Mnemonic	Category/Type							
6907567	Respiratory Virus Mini-Panel by PCR	Mini Vresp PCR	Careset	Careset Component Information						
CPT: 87631	Per the CDC, clinicians are strongly encouraged to test for other causes of respiratory illness (e.g., influenza, RSV, etc) prior to ordering the SARS-CoV-2 PCR test Mini-Flocked Swab in Universal Transport Media (UTM/VTM) or 3mL (1mL) Nasal wash or Bronchial lavage in Sterile Screwtop Container The preferred specimen is mini-Flocked Swab in Universal Transport Media (UTM) (Comes as a kit: RML Supply# 50775), BD Viral Transport Media (VTM) or M4. Keep swabs refrigerated, Stability 7 days refrigerated. Also acceptable 3mL(1mL) BAL or NP/Nasal/Tracheal Aspirate Sterile Screwtop tube Refrigerated. Unacceptable Spcimens: Eswabs, Caclium Alginate or Wooden Shaft Swabs.			Order Code	LOINC	Result Code	Mnemonic	Test Name	Reference Range UOM	Result Type
				5565570			FLU A PCR	Influenza A by PCR		
					34487-9	5565570	Flu A PCR	Influenza A by PCR	N/A	Alpha
				5565580			FLU B PCR	Influenza B by PCR		
					40982-1	5565580	Flu B PCR	Influenza B by PCR	N/A	Alpha
				5565560			RSV PCR	RSV (Respiratory Syncytial Virus) PCR		
					40988-8	5565560	RSV PCR	RSV by PCR	N/A	Alpha
5575875			METAPNEUMO	Metapneumovirus PCR						
	38917-1	5575875	Metapneumo PCR	Metapneumovirus PCR	N/A	Alpha				