

NEW TESTS - Please update your EMR catalog with those appropriate to your practice

NEW TEST Available 12/17/2012									
Order Test Code	Test Name	Mnemonic	Category/Type	Result Code	Result Name	Mnemonic	UOM	Reference Range	Result
3602272	CORTISOL FREE URINE 24H/RANDOM	CORT FR UR	Group Parent Non-Orderable Child	3602222	CORTISOL FREE URINE - UG/24H	XCORT FR U	UG/24H	Adult 24-hr 4.0-50.0 Pediatric 24-hr 1-4 years: 0.9-8.2 5-9 years: 1.0-30.0 10-13 years: 1.0-45.0 3-8 Years: 0.11-0.68 9-12 Years: 0.17-1.41 13-17 Years: 0.29-1.87 Adults: 0.63-2.50	Numeric
			Non-Orderable Child	4504002	CREATININE, 24HR URINE	XCREA/CORT	G/24HR		Numeric
<p>CPT Code: 82530 List Price: \$57.34</p> <p>Specimen Collection: Volume: 2mL (0.5mL) or 24 hr urine container Type: Urine Container: 24 hour Urine, Monovette, Sterile screw top container Transport Environment: Refrigerated Instructions: Adequate refrigeration is the most important aspect of specimen preservation. Collect 24 hour or random Urine. Send refrigerated 24 hr Urine Container or 2mL(0.5mL) aliquot from 24 hour urine collection (re-suspend urine before aliquotting). Keep urine refrigerated during & after collection. Record total urine volume and duration on specimen container and on requisition. STABILITY: Refrigerated=7days, Frozen=30days</p> <p style="text-align: right;">** Possible alpha responses are: See Below, SEE FTN, TNP</p>									

NEW TEST Available 01/14/2013									
Order Test Code	Test Name	Mnemonic	Category/Type	Result Code	Result Name	Mnemonic	UOM	Reference Range	Result
2010900	T3 REVERSE	REVERSE T3	Single Detail	2010900	T3 REVERSE	REVERSE T3	NG/DL	11-32	Numeric
<p>CPT Code: 84482 List Price: \$236.30</p> <p>Specimen Collection: Volume: 1.0mL (0.4mL) Type: Serum Container: Tiger Top Clot Transport Environment: Refrigerated Stability: Refrigerated-5days, Frozen-1year</p> <p style="text-align: center;">****THIS IS A RE-ACTIVATION OF AN EXISTING TEST. IT MAY ALREADY BE IN YOUR SYSTEM ****</p>									

NEW TEST Available 12/17/2012									
Order Test Code	Test Name	Mnemonic	Category/Type	Result Code	Result Name	Mnemonic	UOM	Reference Range	Result
3630327	CHROMATIN (NUCLEOSOMAL) AB	CHROMATIN	Single Detail	3630327	CHROMATIN (NUCLEOSOMAL) AB	CHROMATIN	AI	< 1.0	Numeric
<p>CPT Code: 86325 List Price: \$84.25</p> <p>Specimen Collection: Volume: 1mL (0.5mL) Type: Serum Container: Red No Gel Tube Transport Environment: Refrigerated Stability: Refrigerated-7days; Frozen-30days Instructions: Allow specimen to clot, centrifuge, aliquot and refrigerate.</p>									

MODIFIED TESTS - Please update your EMR catalog with those appropriate to your practice

MODIFICATION Effective Now		
Orderable Test Code	Result Name	Mnemonic
3657650	FELBAMATE (FELBATOL)	FELBAMATE
<p>Update Trough Reference Range: 30-50 mcg/mL</p> <p>Possible Results: Numeric - 0.0-999.9 Alpha - SEE FTN, See Below, TNP</p>		

MODIFICATION Effective Now		
Order Test Code	Test Name	Mnemonic
3618200	TAMBOCOR	TAMBOCOR
<p>Update Specimen Collection: Volume: 2mL(0.5mL) Type: Serum Container: Red No Gel Tube Transport: Refrigerated Stability: Refrigerated-7days, Frozen-30days Instructions: Do not use gel barrier tubes. Optimum time to collect sample: 1 hour before next dose.</p>		

MODIFIED TESTS - Please update your EMR catalog with those appropriate to your practice

MODIFICATION		Effective 12/17/2012	
Order Test Code	Test Name	Mnemonic	
5508985	FEBRIL AGG (SALMONELLA T. Abs)	FEBRIL AGG	
			<p>Update Reference Range: Negative</p> <p>Possible Results: NEG, POS, SEE FTN, TNP</p> <p>Update Specimen Collection: Volume: 1mL(0.1mL) Type: Serum Container: Serum Separator or Red No Gel Tube Transport: Refrigerated Stability: Refrigerated-7days, Frozen-30days Instructions: Order must be approved by RML doctor.</p>

MODIFICATION		Effective 01/14/2013	
Order Test Code	Test Name	Mnemonic	
2004240	MYOGLOBIN SERUM	MYOGLOBIN	
			<p>Update Specimen Collection: Volume: 1mL(0.5mL) Type: Serum Container: Tiger Clot tube Transport: Refrigerated Stability: Refrigerated-7days, Frozen-35days Instructions: Plasma is unacceptable.</p>

MODIFICATION		Effective 12/17/2012	
Order Test Code	Test Name	Mnemonic	
5574700	ANTI PARVOVIRUS B-19 IGG/IGM	PARVO B19	
			<p>Update Reference Ranges for non-orderable result details 5574600 and 557465: <0.9 Negative 0.9-1.1 Equivocal >1.1 Positive</p> <p>Remove non-orderable result detail: 5574750</p> <p>Results: Numeric = 0.0-99.9 or Alpha= See Below, SEE FTN or TNP</p> <p>Update Specimen Collection: Volume: 1mL(0.2mL) Type: Serum Container: Serum Separator or Red No Gel Tube Transport: Refrigerated Stability: Refrigerated-7days, Frozen-30days Instructions: Please mark specimen plainly as "acute" or "convalescent".</p>

MODIFICATION		Effective	
Order Test Code	Test Name	Mnemonic	
3653300	AMIODARONE LEVEL	AMIODARON	
			<p>EFFECTIVE NOW</p> <p>Update Specimen Collection: Volume: 3mL(1mL) Type: Serum Container: Red Top No Gel Tube Transport: Refrigerated Stability: Refrigerated-4days, Frozen-30days Instructions: Allow to clot, centrifuge, separate from cells immediately and refrigerate.</p> <p>EFFECTIVE 01/14/2013</p> <p>Update CPT to: 80299 x2</p> <p>Update Reference Range: 1.5-2.5 mcg/mL</p>

INACTIVE TESTS - Please update your EMR catalog with those appropriate to your practice

INACTIVE		Effective 12/17/2012		
Order Test Code	Test Name	Mnemonic	Replacement Code	Replacement Name
3602275	CORTISOL FREE URINE 24hr/RANDM	CORT FR U	3602272	CORTISOL FREE URINE 24H/RANDOM

INACTIVE		Effective 01/14/2013		
Order Test Code	Test Name	Mnemonic	Replacement Code	Replacement Name
2010915	REVERSE T3	REVERSE T3	2010900	T3 REVERSE