



All specimens must be labeled with a first and last name and a 2nd identifier

RML SPECIMEN COLLECTION TRANSPORT CONTAINERS


Revised 5/9/2018




Aptima Multi-test swab- Vaginal Swab-Chlamydia, Gonorrhoeae, Bacterial Vaginosis, Candida (5 Species), Herpes, Trichomonas



Aptima Unisex- Endocervical Swab-Chlamydia, Gonorrhoeae, Trichomonas



Aptima Urine-Male or Female-Chlamydia, Gonorrhoeae, Trichomonas



Copan ESwab

- *Aerobic Cultures
- *Aerobic Wound Cultures
- *Anaerobic Wound Cultures
- *Ear, Nose & Throat Cultures
- *Fungal Cultures
- *Genital/Vaginal Cultures (GBS)
- *Rectal Cultures
- Strep Cultures
- Rapid Strep A
- Group B Strep



Green Copan ESwab


- *Eye Cultures
- *Urethral Cultures
- *Nasopharyngeal Cultures
- *MRSA SCREEN



Histo Container

Thinprep- PAP


- *HPV
- *Chlamydia
- *Gonorrhoeae
- *Trichomonas



Virus Culture & Chlamydia Culture (3 mL)


- *Mycoplasma
- *Ureaplasma
- *Herpes Culture
- *Rapid Influenza (Flu A/FluB)
- *Peds RSV (1mL)
- *VRESP PCR
- *Bordetella pertussis/parapertussis
- * Plastic tube containing pink virus media and a swab
- * Break off swab inside tube and replace cap

UTM * **Transport on ice or refrigerate**



Wet Prep Saline Media

- * Wet prep for yeast and Trichomonas
- * Plastic tube containing clear saline and a swab
- * Break swab inside container & replace cap



BD Affirm Ambient Temperature Transport Swabs (ATTS)

- *Gardnerella
- *Candida (7 species)
- *Trichomonas


VAG MOLEC



Bactec Bottles for Blood Cultures

- *Blue-Aerobic
- *Yellow-Anerobic
- *Pink-Peds
- *Red-AFB/Fungus
- *Gray-Antibiotics

(Do not label over Barcode on bottle.)




Sterile Cup

Body Fluids

- *Respiratory Cultures
- *Urine Cultures
- *Stool for Polys or leukocytes
- *Fungal Cultures (skin scrapings, hair, nails)
- *AFB Cultures

Tissues

- *Skin Scrapings for KOH
- *C Diff AN (if performed in 24hr)



ETM

- *C. Diff AN (up to 5 day)
- *Stool Culture-ETM




Pin Worm Prep



PARAPAK

- *Stool for Ova & Parasites
- *Stool for Polys or leukocytes



ICT OCCULT VIAL

Used when Physician orders ICT OCCULT