

NEW REPORT - Please update your EMR catalog with those appropriate to your practice

NEW REPORT Effective : 02/04/2019									
Test Code	REPORT Name	Mnemonic	Category/Type	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
	FLOW CYTOMETRY FINAL REPORT	NA	Report	7120723	Flow Cytometry Final Report	Flow Cytometry Final Report	NA	NA	
Note: This is a new AP Report Result									

NEW TESTS - Please update your EMR catalog with those appropriate to your practice

New Test Effective : 02/18/2019									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6906783	Breast & Ovarian Hereditary Cancer Panel, Sequencing & Deletion/Duplication, 20	Breast/Ovarian Cancer Pnl	Group	35474-6	6906785	Brst/OvarCAIntp	Brst/OvarCAIntp	By Report	Freetext
				31208-2	6906787	Brst/OvarCA Spc	Brst/OvarCA Spc	By Report	Freetext
CPT: 81432	Specimen Prep: 3mL (1mL) Whole Blood. Send Refrigerated.								
81433	Collect: Lavendar or Yellow (ACD solution A or B) Submit with order: Pateint history form Stability: Room Temperature - 72 hours; Refrigerated - 1 week; Frozen - Unacceptable Method: Massively Parallel Sequencing/Exonic Oligonucleotide-based CGH Microarry Set-up: Varies TAT Varies: 6-8 weeks set-up								

New Test Effective : 02/18/2019

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6906815	Breast & Ovarian Hearedity Cancer Syndrome (BRCA1/BRCA2) Sequencing & Deletion/Duplication	BRCA1/BCRA2 Seq/Del/Dup	Group	59041-4	6906817	BRCA1/2 Intrp	BRCA1/2 Intrp	By Report	Freetext
				31208-2	6906819	BRCA1/2 Spec	BRCA1/2 Spec	By Report	Freetext
CPT: 81162	Specimen Prep: 3mL (1mL) Whole Blood. Send Refrigerated. Collect: Lavendar or Yellow (ACD solution A or B) Submit with order: Pateint history form Stability: Room Temperature - 72 hours; Refrigerated - 1 week; Frozen - Unacceptable Method: Polymerase Chain Reaction/Sequencing/Multiplex Ligation-dependent Probe Amplification Set-up: Sun-Sat TAT Varies: 14-21 day after set-up								

New Test Effective : 02/18/2019 Updated Information

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6906789	Breast & Ovarian Hearedity Cancer Syndrome (BRCA1/BRCA2) Sequencing	BRCA1/BCRA2 Sequencing	Group	59041-1	6906791	BRCA1/2 Seq Intrp	BRCA1/2 Intrp	By Report	Freetext
				31208-2	6906793	BRCA1/2 Seq Spec	BRCA1/2 Spec	By Report	Freetext
CPT: 81163	Specimen Prep: 3mL (1mL) Whole Blood. Send Refrigerated. Collect: Lavendar or Yellow (ACD solution A or B) Submit with order: Pateint history form Stability: Room Temperature - 72 hours; Refrigerated - 1 week; Frozen - Unacceptable Method: PCR/Sequencing Set-up: Sun-Sat TAT Varies: 14-21 day after set-up								
<i>Updated to correct CPT code from 81211 to 81163</i>									

MODIFICATION TESTS - Please update your EMR catalog with those appropriate to your practice

Modify Test Effective : 2/4/19 - Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4003230	Naproxen	Naproxen	Detail	N/A	4003230	Naproxen	Naproxen		
1. Minimum volume increase from 0.21mL to 0.25mL 2. Stability updated: Ambient/refrigerated/frozen - 15 days 3. Testing Schedule Changes: Current Set-Up: Tues & Thurs; TAT: 3 Days New Set-up: Daily; TAT: 7 days 4. Methodology change from HPLC to HPLC/Tandem Mass Spec. No other changes									

Modify Test Effective : 2/11/19 - Please update your EMR catalog with those appropriate to your practice -									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
3001875	Thallium 24 Hour Urine	THALIU 24U	Detail	NA	3001875	THALIU 24U	Thallium 24 Hour Urine		
<p>1. TAT updated: From 1 day after set-up to 2-5 days after set-up 2. Testing schedule will not change. No other changes</p>									

Modify Test Effective : IMMEDIATELY 12/10/18 - Please update your EMR catalog with those appropriate to your practice -																																																																																																		
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type																																																																																									
3603055	Deoxycorticosterone Level	DEOXYSTER LVL	Detail	NA	3603055	Deoxyster Lvl	Deoxycorticosterone Level	See Below																																																																																										
<p>1. Unacceptable specimen types: Grossly lipemic, grossly icteric 2. Method Changing from LC/MS/MS to Chromatography/Mass Spectrometry 3. Reference Range Changes:</p>																																																																																																		
<table border="1"> <thead> <tr> <th colspan="2">Current Ref Range</th> </tr> </thead> <tbody> <tr> <td>Male ≥18 Years</td> <td>≤15 ng/dL</td> </tr> <tr> <td>Female-Phase of Menstrual Cycle</td> <td></td> </tr> <tr> <td> Mid Follicular</td> <td>≤18 ng/dL</td> </tr> <tr> <td> Surge</td> <td>≤23 ng/dL</td> </tr> <tr> <td> Mid Luteal</td> <td>≤19 ng/dL</td> </tr> <tr> <td>Pediatric 16-17 Years</td> <td>≤35 ng/dL</td> </tr> </tbody> </table>				Current Ref Range		Male ≥18 Years	≤15 ng/dL	Female-Phase of Menstrual Cycle		Mid Follicular	≤18 ng/dL	Surge	≤23 ng/dL	Mid Luteal	≤19 ng/dL	Pediatric 16-17 Years	≤35 ng/dL	<table border="1"> <thead> <tr> <th colspan="3">New Ref Ranges</th> </tr> <tr> <th>Age</th> <th>Male</th> <th>Female</th> </tr> </thead> <tbody> <tr> <td><30 days</td> <td>No Range Established</td> <td>No Range Established</td> </tr> <tr> <td>1-11 months</td> <td>< or = 69 ng/dL</td> <td>< or = 93 ng/dL</td> </tr> <tr> <td>1 year</td> <td>< or = 66 ng/dL</td> <td>< or = 100 ng/dL</td> </tr> <tr> <td>2 years</td> <td>< or = 64 ng/dL</td> <td>< or = 100 ng/dL</td> </tr> <tr> <td>3 years</td> <td>< or = 63 ng/dL</td> <td>< or = 96 ng/dL</td> </tr> <tr> <td>4 years</td> <td>< or = 63 ng/dL</td> <td>< or = 89 ng/dL</td> </tr> <tr> <td>5 years</td> <td>< or = 62 ng/dL</td> <td>< or = 78 ng/dL</td> </tr> <tr> <td>6 years</td> <td>< or = 58 ng/dL</td> <td>< or = 58 ng/dL</td> </tr> <tr> <td>7 years</td> <td>< or = 51 ng/dL</td> <td>< or = 53 ng/dL</td> </tr> <tr> <td>8 years</td> <td>< or = 44 ng/dL</td> <td>< or = 45 ng/dL</td> </tr> <tr> <td>9 years</td> <td>< or = 38 ng/dL</td> <td>< or = 40 ng/dL</td> </tr> <tr> <td>10 years</td> <td>< or = 31 ng/dL</td> <td>< or = 35 ng/dL</td> </tr> <tr> <td>11 years</td> <td>< or = 26 ng/dL</td> <td>< or = 30 ng/dL</td> </tr> <tr> <td>12 years</td> <td>< or = 22 ng/dL</td> <td>< or = 25 ng/dL</td> </tr> <tr> <td>13 years</td> <td>< or = 20 ng/dL</td> <td>< or = 21 ng/dL</td> </tr> <tr> <td>14 years</td> <td>< or = 20 ng/dL</td> <td>< or = 19 ng/dL</td> </tr> <tr> <td>15 years</td> <td>< or = 20 ng/dL</td> <td>< or = 18 ng/dL</td> </tr> <tr> <td>16 years</td> <td>< or = 20 ng/dL</td> <td>< or = 18 ng/dL</td> </tr> <tr> <td>17 years</td> <td>< or = 19 ng/dL</td> <td>< or = 18 ng/dL</td> </tr> <tr> <td>Adult</td> <td>< or = 18 ng/dL</td> <td></td> </tr> <tr> <td> Adult - Pre-menopausal Mid Follicular</td> <td>< or = 18 ng/dL</td> <td></td> </tr> <tr> <td> Adult - Pre-menopausal Surge</td> <td>< or = 23 ng/dL</td> <td></td> </tr> <tr> <td> Adult - Pre-menopausal Mid Luteal</td> <td>< or = 19 ng/dL</td> <td></td> </tr> </tbody> </table>						New Ref Ranges			Age	Male	Female	<30 days	No Range Established	No Range Established	1-11 months	< or = 69 ng/dL	< or = 93 ng/dL	1 year	< or = 66 ng/dL	< or = 100 ng/dL	2 years	< or = 64 ng/dL	< or = 100 ng/dL	3 years	< or = 63 ng/dL	< or = 96 ng/dL	4 years	< or = 63 ng/dL	< or = 89 ng/dL	5 years	< or = 62 ng/dL	< or = 78 ng/dL	6 years	< or = 58 ng/dL	< or = 58 ng/dL	7 years	< or = 51 ng/dL	< or = 53 ng/dL	8 years	< or = 44 ng/dL	< or = 45 ng/dL	9 years	< or = 38 ng/dL	< or = 40 ng/dL	10 years	< or = 31 ng/dL	< or = 35 ng/dL	11 years	< or = 26 ng/dL	< or = 30 ng/dL	12 years	< or = 22 ng/dL	< or = 25 ng/dL	13 years	< or = 20 ng/dL	< or = 21 ng/dL	14 years	< or = 20 ng/dL	< or = 19 ng/dL	15 years	< or = 20 ng/dL	< or = 18 ng/dL	16 years	< or = 20 ng/dL	< or = 18 ng/dL	17 years	< or = 19 ng/dL	< or = 18 ng/dL	Adult	< or = 18 ng/dL		Adult - Pre-menopausal Mid Follicular	< or = 18 ng/dL		Adult - Pre-menopausal Surge	< or = 23 ng/dL		Adult - Pre-menopausal Mid Luteal	< or = 19 ng/dL	
Current Ref Range																																																																																																		
Male ≥18 Years	≤15 ng/dL																																																																																																	
Female-Phase of Menstrual Cycle																																																																																																		
Mid Follicular	≤18 ng/dL																																																																																																	
Surge	≤23 ng/dL																																																																																																	
Mid Luteal	≤19 ng/dL																																																																																																	
Pediatric 16-17 Years	≤35 ng/dL																																																																																																	
New Ref Ranges																																																																																																		
Age	Male	Female																																																																																																
<30 days	No Range Established	No Range Established																																																																																																
1-11 months	< or = 69 ng/dL	< or = 93 ng/dL																																																																																																
1 year	< or = 66 ng/dL	< or = 100 ng/dL																																																																																																
2 years	< or = 64 ng/dL	< or = 100 ng/dL																																																																																																
3 years	< or = 63 ng/dL	< or = 96 ng/dL																																																																																																
4 years	< or = 63 ng/dL	< or = 89 ng/dL																																																																																																
5 years	< or = 62 ng/dL	< or = 78 ng/dL																																																																																																
6 years	< or = 58 ng/dL	< or = 58 ng/dL																																																																																																
7 years	< or = 51 ng/dL	< or = 53 ng/dL																																																																																																
8 years	< or = 44 ng/dL	< or = 45 ng/dL																																																																																																
9 years	< or = 38 ng/dL	< or = 40 ng/dL																																																																																																
10 years	< or = 31 ng/dL	< or = 35 ng/dL																																																																																																
11 years	< or = 26 ng/dL	< or = 30 ng/dL																																																																																																
12 years	< or = 22 ng/dL	< or = 25 ng/dL																																																																																																
13 years	< or = 20 ng/dL	< or = 21 ng/dL																																																																																																
14 years	< or = 20 ng/dL	< or = 19 ng/dL																																																																																																
15 years	< or = 20 ng/dL	< or = 18 ng/dL																																																																																																
16 years	< or = 20 ng/dL	< or = 18 ng/dL																																																																																																
17 years	< or = 19 ng/dL	< or = 18 ng/dL																																																																																																
Adult	< or = 18 ng/dL																																																																																																	
Adult - Pre-menopausal Mid Follicular	< or = 18 ng/dL																																																																																																	
Adult - Pre-menopausal Surge	< or = 23 ng/dL																																																																																																	
Adult - Pre-menopausal Mid Luteal	< or = 19 ng/dL																																																																																																	
<p>No other changes</p>																																																																																																		

Modify Test Schedule Effective : 2/25/19 - Please update your EMR catalog with those appropriate to your practice -									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
5569625	Candida Albicans Antibody Panel	CANDIDA ABS	Group	27391-2	5569925	C Albicans IgM	Candida Albicans IgM		
				07172-0	5569675	C Albicans IgG	Candida Albicans IgG		
				27417-5	5570025	C Albicans IgA	Candida Albicans IgA		
<p>1. Test Set-Up Changes: Current: Setup-Tues; TAT 1-8 days NEW: Setup-Mon, Wed, Sat; TAT 1-3 days after setup. No other changes</p>									

INACTIVATED TESTS - Please update your EMR catalog with those appropriate to your practice

INACTIVE Effective : 01/01/2018									
Inactive Test Code	Inactive Test Name	Inactive Mnemonic	Replacement Information						
			Test Code	Mnemonic	Test Name				
4307075	Drug Screen, Breast Milk	DRUG BRMLK	4315450	Mlk Drg Sc	Drug Screen Breast Milk				
			4315325	B Mlk Amphet	Amphetamines Breast Milk				
			4315350	B Mlk Barb	Barbiturates Breast Milk				
			4315375	B Mlk Benzo	Benzodiazapines Breast Milk				
			4315400	B Mlk Cocaine	Cocaine Breast Milk				
			4315425	B Mlk Methad	Methadone Breast Milk				
			4315475	B Mlk PCP	Phencyclidine (PCP) Breast Milk				
			4315500	B Mlk Propox	Propoxyphene Breast Milk				
			4315525	B Mlk THC	Tetrahydrocannabinol (THC) Breast Milk				
			4315550	B Mlk Opiate	Opiates Breast Milk				
<p>CPT: This test is being discontinued due to low utilization</p>									