







Bone Marrow Collection Quick Reference Guide

Every Bone Marrow case should have the following specimens collected:

Whole Blood	Bone Marrow Aspirate	Bone Marrow Clot	Bone Marrow Core Biopsy	Other
<p>2 EDTA Tubes</p> 	<p>1 EDTA Tube + 1 Na Citrate + 1 Na Heparin</p>  <p>(2 mL min per tube)</p>	<p>Formalin Container</p> 	<p>Formalin Container</p>  <p>Dry Tap Collections</p> <p>Prepare 8-10 touch preps, section the biopsy into 3 pieces and place 1 section into each container:</p> <ul style="list-style-type: none"> 1 EDTA tube w/RPMI 1 Na Heparin tube w/ RPMI 1 Formalin Container 	<p>Paperwork</p> <ol style="list-style-type: none"> General Requisition Hematopathology Requisition Bone Marrow Collection Form <p>SPS/Isolator Tube (Optional for Microbiology)</p> 
<p>Label all specimens with patient name, MRN, DOB, collection date and time AND the specimen type (i.e. blood, aspirate, clot, biopsy)</p>				