



MEMORANDUM

To: Regional Medical Laboratory Clients
 From: Regional Medical Laboratory, Inc.
 Sonja Matthews, MT, (ASCP) SC, Director of Lab Operations
 Re: **Blood Collection Tube Conversion**
 Date: **November 18, 2019**

Please be advised that Regional Medical Laboratory, Inc. will be converting blood collection tubes from BD to Greiner starting November 18, 2019. Online and faxed supply orders through the RML Materials Management department will automatically be converted for clients. RML will continue to accept BD tubes until current supplies are depleted. RML has completed verification studies per manufacturer guidelines and there are no operational or testing impacts anticipated from this change.

BD Part Number	BD Description	Greiner Part Number	Tube Description
363083	2.7 mL light blue	454334	13x75x3 mL 9NC Coagulation
367815	6 mL red no gel	456089	13x100x5 mL 6mL Serum Clot
367856	3 mL short draw lavender	454246	13x75x3 mL K2EDTA K2 Lavender
367861	4 mL lavender	454209	13x75x4 mL K2EDTA K2 Lavender
367878	6 mL Na heparin	456028	13x100x6 mL NH Sodium Heparin
367886	6 mL lithium hep no gel	456088	13x100x6 mL LH Lithium Heparin
367899	6 mL pink	456279	13x100x6 mL K2 EDTA K2 Pink
367922	4 mL gray Na fluoride	454297	13x75x4 mL Vacuette Glycolytic
367962	4.5 mL lithium hep gel	456087P	13x100x5 mL LH Lithium Heparin
367986	5 mL gold gel	456018P	13x100x5 mL Serum Clot Activator
368380	6 mL royal no additive	456275	13x100x6 mL NH Sodium Heparin
368381	6 mL royal trace metal EDTA	456275	13x100x6 mL NH Sodium Heparin

Please contact your RML Account Representative, Client Services or Sonja Matthews if you have questions regarding this conversion. 918-744-2500 or 1-800-722-8077