

NEW TESTS - Please update your EMR catalog with those appropriate to your practice

New Test Effective : 04/28/2020

| Test Code | Test Name | Mnemonic | Category/Type | LOINC | Result Code | Mnemonic | Result Name | Reference Range UOM | Result Type |
|-----------|----------------------------|----------------|---------------|---------|-------------|----------------|----------------------------|---------------------|-------------|
| 6901550 | SARS-CoV-2 IgG Qualitative | SARS-CoV-2 IgG | Detail | 94507-1 | 6901550 | SARS-CoV-2 IgG | SARS-CoV-2 IgG Qualitative | Negative | Alpha |

CPT: 86769

aka: COVID-19 IgG Ab; Coronavirus 2019 IgG Antibody

Preferred specimen: 2 mL (0.5 mL) Serum in Clot Activator SST or EDTA Plasma Aliquot
 Storage/Transport: Refrigerated
 Stability: Ambient 2 days; Refrigerated 7 days; Freeze if testing will be delayed more than 7 days
 Testing Schedule: Mon-Fri / Expected TAT: 1-2 Days
 Methodology: Chemiluminescent Microparticle Immunoassay
 Clinical Use: Detects IgG antibodies against the 2019 novel coronavirus (SARS-CoV-2/COVID-19) to evaluate exposure. This test should not be used for screening of donated blood