

NEW TESTS - Please update your EMR catalog with those appropriate to your practice

New Test Effective : 11/16/2020									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6907511	BRCA1/2 NGS Sequencing and Del/Dup	BRCA1/2 NGS S/D/D	Group						
CPT: 81162	Collect: Lavender or Yellow ACD Solution A or B Specimen: 3mL (2mL) EDTA whole blood Refrigerated. DO NOT SPIN OR FREEZE! Unacceptable: Serum or plasma; grossly hemolyzed or frozen specimens Stability: RT: 72 hours; RF: 1 week; FZ: Unacceptable Method: MPS/MLPA Set-up: Varies; TAT: 15-17 days after set-up			31208-2	6907513	BRCA1/2 Spec	BRCA1/2 Specimen	By report	Freertext
				94191-4	6907515	BRCA1/2 Interp	BRCA1/2 Interpretation	By report	Freertext

New Test Effective : 11/16/2020									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6907491	Coccidioides Ab Reflexive Panel	Coccid Ab Rflx Pnl	Group						
CPT: 86635	Collect: Serum Separator Tube (SST). Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Parallel testing is preferred and convalescent specimens must be received within 30 days from receipt of the acute specimens. Storage/Transport Temperature: Refrigerated. Remarks: Please mark specimens plainly as "acute" or "convalescent." Unacceptable Conditions: Contaminated, hemolyzed, or severely lipemic specimens. Stability (collection to initiation of testing): Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year (avoid repeated freeze/thaw cycles) Set-up: Sun - Sat; TAT: 1-6 days after set-up <i>See Reflex Testing below</i>			27009-0	6907493	Coccid IgG	Coccidioides IgG	XX.XX	Free text
				90448-2	6907495	Coccid IgM	Coccidioides IgM	XX.XX	Free text

Reflex Test Name	Reflex Mnemonic	Category/Type	Order Code	LOINC	Result Code	Mnemonic	Test Name	Reference Range UOM	Result Type
.Coccidioides Titer by CF	.Coccid Titer	Detail							
CPT: 86635 NOT ORDERABLE ALONE - ONLY AS REFLEX This is a reflex test to Coccidioides Ab Reflexive Panel (6907491) REFLEX CRITERIA: Negative result for Coccidioides Ab Reflexive Panel (6907491)			6907497			.Coccid Titer	.Coccidioides Titer by CF		
			33380-7	6907497	Coccid Titer	Coccidioides Titer by CF	XX.XX	Free text	

Reflex Test Name	Reflex Mnemonic	Category/Type	Order Code	LOINC	Result Code	Mnemonic	Test Name	Reference Range UOM	Result Type
.Coccidioides Ab by ID	.Coccid Ab ID	Detail							
CPT: 86635 NOT ORDERABLE ALONE - ONLY AS REFLEX This is a reflex test to Coccidioides Ab Reflexive Panel (6907491) REFLEX CRITERIA: Negative result for Coccidioides Ab Reflexive Panel (6907491)			6907499			.Coccid Ab ID	.Coccidioides Ab by ID		
			5095-5	6907499	.Coccid Ab ID	Coccidioides Ab by ID	XX.XX	Free Text	

New Test Effective : 11/16/2020									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6908061	Babesia duncani (WA1) IgG Ab	B duncani IgG	Detail						
CPT: 86753	Specimen: 0.5mL (0.2mL) serum Temperature: Refrigerated Method: IFA Stability: RT-7 days; RF-14 days; FZ-30 days Set-up: Tue-Sat; TAT: 1-3 days after set-up			41415-1	6908061	B duncani IgG	Babesia duncani (WA1) IgG Ab	XXXXXX titer	Freertext

MODIFICATION TESTS - Please update your EMR catalog with those appropriate to your practice

Modify Test Effective : Immediately - Please update your EMR catalog with those appropriate to your practice -									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4400350	Cortisol Serum Total	TOT CORTIS	Detail						
CPT: 82533	To Replace the Discontinued Cortisol AM and Cortisol PM tests the Cortisol Total will be reported with a footnote containing the reference ranges for AM and PM collections. For specific timed collections: The AM cortisol level should be collected within 1-2 hours after the normal waking time of the patient. The PM cortisol level should be collected within 8-14 hours after the normal waking time of the patient. Stability: Ambient 1 day. Refrigerated 14 days. Frozen 30 days. Methodology: Chemiluminescence Assay Testing Schedule Sun - Fri / TAT 1-3 Days			43215-3	4501800	Cortisol	Cortisol	by report ug/dL	Numeric XXXX.X

Modify Test		Effective : Immediately		- Please update your EMR catalog with those appropriate to your practice -							
Test Code	Test Name	Mnemonic	Category/Type	Careset Component Information							
6908087	ENT Allergy Panel 1 (Large panel)	ENT PNL 1	Careset	Order Code	LOINC	Result Code	Mnemonic	Test Name	Reference Range	UOM	Result Type
3 careset components added 3 careset components removed No other test changes				5606795			Elm Tree	Elm Tree Allergen			
					15697-6	5506785	Elm Tree CL	Elm Tree Allergen	See below	See below	
					6109-3	5606775	Elm Tree CT	Elm Tree Count	See below	See below	
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below	
					5613000		Silver Birch Tree	Common Silver Birch Tree Allergen			
					15283-5	5612075	Silver Birch CL	Common Silver Birch Tree Allergen	See below	See below	
					15579-6	5612050	Silver Birch CT	Common Silver Birch Tree Count	See below	See below	
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below	
					5611750		Box Elder	Box Elder (Maple) Tree Allergen			
					15585-3	5611725	Box Elder CL	Box Elder Tree Allergen	See below	See below	
					7155-5	5611700	Box Elder CT	Box Elder Tree Count	See below	See below	
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below	
					5607475		Cottonwood Tree	Cottonwood Tree Allergen			
					15659-6	5607450	Cottonwood CL	Cottonwood Tree Allergen	See below	See below	
					6090-5	5607425	Cottonwood CT	Cottonwood Tree Count	See below	See below	
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below	
					5518175		Mulberry Tree Red	Mulberry Tree Red Allergen			
					21400-7	5518150	Mulberry Red CL	Mulberry Tree Red Allergen	See below	See below	
					7514-3	5528650	Mulberry Red CT	Mulberry Tree Red Count	See below	See below	
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below	
					5528010		Mulberry Tree White	Mulberry Tree White Allergen			
					N/A	5528005	Mulberry Wht CT	Mulberry Tree White Allergen	See below	See below	
					N/A	5528015	Mulberry Wht CL	Mulberry Tree White Count	See below	See below	
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below	
					5611850		Mountain Juniper	Mountain Juniper Tree Allergen			
					15615-8	5611825	Mtn Juniper CL	Mountain Juniper Tree Allergen	See below	See below	
					6178-8	5611800	Mtn Juniper CT	Mountain Juniper Tree Count	See below	See below	
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below	
					5604925		Oak Tree	Oak Tree Allergen			
					N/A	5590800	Oak Tree CL	Oak Tree Allergen	See below	See below	
					50030-6	5528800	Oak Tree CT	Oak Tree Count	See below	See below	
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below	
					5606750		Pecan, Hickory	Pecan, Hickory Tree Allergen			
					15774-3	5606725	Pecan Tr CL	Pecan, Hickory Tree Allergen	See below	See below	
					6209-1	5606700	Pecan Tr CT	Pecan, Hickory Tree Count	See below	See below	
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below	
					5621250		Privet Tree	Privet Tree Allergen			
					15960-8	5619050	Privet Tree CL	Privet Tree Allergen	See below	See below	
					7632-3	5619025	Privet Tree CT	Privet Tree Count	See below	See below	
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below	
		5613519		White Pine	White Pine Tree IgE						
		6282-8	5613519	White Pine Tr IgE	White Pine Tr IgE	See below	See below				
		5613225		Mesquite Tree	Mesquite Tree Allergen						
		15843-6	5613200	Mesquite Tr CL	Mesquite Tree Allergen	See below	See below				
		6173-9	5613175	Mesquite Tr CT	Mesquite Tree Count	See below	See below				
		8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below				
		5606475		White Ash	White Ash Tree Allergen						
		15546-5	5606450	Wht Ash Tr CL	White Ash Tree Allergen	See below	See below				
		6278-6	5606425	Wht Ash Tr CT	White Ash Tree Count	See below	See below				
		8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below				
		5613499		Red Cedar	Cedar Red Tree IgE						
		7195-1	5613499	Red Cedar Tr IgE	Red Cedar Tr IgE	See below	See below				
		5606325		Bermuda Grass	Bermuda Grass Allergen						
		15738-8	5606300	Bermuda CL	Bermuda Grass Allergen	See below	See below				
		6041-8	5606275	Bermuda CT	Bermuda Grass Count	See below	See below				
		8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below				
		5605965		Johnson Grass	Johnson Grass Allergen						
		15743-8	5605950	Johnson Gr CL	Johnson Grass Allergen	See below	See below				
		6152-3	5605925	Johnson Gr CT	Johnson Grass Count	See below	See below				
		8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below				
		5611950		Timothy Grass	Timothy Grass Allergen						
		16054-9	5611925	Timothy Gr CL	Timothy Grass Allergen	See below	See below				
		6265-3	5611900	Timothy Gr CT	Timothy Grass Count	See below	See below				
		8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below				

6908087 ENT Allergy Panel 1 (Large panel) - - Continued on next page

6908087 ENT Allergy Panel 1 (Large panel) - - -Continued from previous page

5617075			Bahia Grass	Bahia Grass Allergen		
	15737-0	5617050	Bahia Grass CL	Bahia Grass Allergen	See below	See below
	6034-3	5617025	Bahia Grass CT	Bahia Grass Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5609225			Cocklebur	Cocklebur Allergen		
	15646-3	5609200	Cocklebur CL	Cocklebur Allergen	See below	See below
	6077-2	5609175	Cocklebur CT	Cocklebur Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5582300			Dandelion	Dandelion Allergen		
	15676-0	5582250	Dandelion CL	Dandelion Allergen	See below	See below
	6097-0	5582200	Dandelion CT	Dandelion Count	See below	See below
5607025			Giant Ragweed	Giant Ragweed Allergen		
	15978-0	5607000	Giant Ragweed CL	Giant Ragweed Allergen	See below	See below
	6124-2	5606975	Giant Ragweed CT	Giant Ragweed Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5621200			Firebush (Kochia)	Firebush (Kochia) Allergen		
	N/A	5621225	Firebush CL	Firebush (Kochia) Allergen	See below	See below
	7491-4	5621175	Firebush CT	Firebush (Kochia) Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5609300			Goosefoot Lamb's Quarte	Lamb's Quarters Goosefoot Allergen		
	15805-5	5609275	Lamb Quarter CL	Lamb's Quarters Goosefoot Allergen	See below	See below
	6156-4	5609250	Lamb Quarter CT	Lamb's Quarters Goosefoot Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5606835			Rough Marshelder	Rough Marshelder Allergen		
	15693-5	5606825	Rough Marsh CL	Rough Marshelder Allergen	See below	See below
	6232-3	5606800	Rough Marsh CT	Rough Marshelder Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5644025			Common Pigweed	Common Pigweed Allergen		
	15936-8	5644075	Pigweed CL	Common Pigweed Allergen	See below	See below
	7604-2	5644050	Pigweed CT	Common Pigweed Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5607625			Plantain (English)	Plantain (English), Ribwort Allergen		
	15952-5	5607600	Plantain Rib CL	Plantain (English), Ribwort Allergen	See below	See below
	6110-1	5607575	Plantain Rib CT	Plantain (English), Ribwort Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5607100			Russian Thistle	Saltwort, Russian Thistle Allergen		
	16052-3	5607072	Russian This CL	Saltwort, Russian Thistle Allergen	See below	See below
	6234-9	5607050	Russian This CT	Saltwort, Russian Thistle Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5609075			Sheep Sorrel	Sheep Sorrel Allergen		
	16016-8	5609050	Sheep Sorrel CL	Sheep Sorrel Allergen	See below	See below
	6244-8	5609025	Sheep Sorrel CT	Sheep Sorrel Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5604625			Common Ragweed	Common Ragweed Allergen		
	15975-6	5604600	Ragweed CL	Common Ragweed Allergen	See below	See below
	6085-5	5604575	Ragweed CT	Common Ragweed Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5650125			Nettle	Nettle Allergen		
	15872-5	5650075	Nettle CL	Nettle Allergen	See below	See below
	6186-1	5650025	Nettle CT	Nettle Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5612025			Mugwort	Mugwort Allergen		
	15863-4	5612000	Mugwort CL	Mugwort Allergen	See below	See below
	6183-8	5611975	Mugwort CT	Mugwort Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5606250			Cat Dander	Cat Dander Allergen		
	15609-1	5606225	Cat Dander CL	Cat Dander Allergen	See below	See below
	6833-8	5606200	Cat Dander CT	Cat Dander Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5613349			Cockroach American	Cockroach American IgE		
	30170-5	5613349	Cockroach Amer IgE	Cockroach American IgE	See below	See below
5606025			German Cockroach	German Cockroach Allergen		
	15647-1	5606000	Gr Roach CL	German Cockroach Allergen	See below	See below
	24139-8	5605975	Gr Roach CT	German Cockroach Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5605800			Dog Dander	Dog Dander Allergen		
	15685-1	5605750	Dog Dander CL	Dog Dander Allergen	See below	See below
	6098-8	5605775	Dog Dander CT	Dog Dander Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below

6908087 ENT Allergy Panel 1 (Large panel) - - - Continued on next page

6908087 ENT Allergy Panel 1 (Large panel) - -Continued from previous page

5604985			Dust Mite (D. farinae)	Dust Mite (D. farinae) Allergen		
	15680-2	5604975	Dust Mite CL	House Dust Mite Allergen	See below	See below
	6095-4	5604950	Dust Mite CT	House Dust Mite Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5607550			House Dust Mite	House Dust Mite Allergen		
	15682-8	5607525	House Dust D Petrony M	House Dust Mite D. Petrony Class	See below	See below
	6096-2	5607500	House Dust D Petrony M	House Dust Mite D. Petrony Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5607400			Feathers	Feathers Allergen		
	31161-3	5607350	Feathers CT	Feathers Count	See below	See below
	24495-4	5607375	Feathers CL	Feathers Allergen	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5606550			Alternaria Alternata	Alternaria Alternata Allergen		
	15530-9	5606525	Alt Ten CL	Alternaria Alternata Allergen	See below	See below
	6020-2	5606500	Alt Ten CT	Alternaria Alternata Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5606400			Aspergillus fumigatus	Aspergillus fumigatus Allergen		
	15549-9	5606375	Asp Fum CL	Aspergillus fumigatus Allergen	See below	See below
	6025-1	5606350	Asp Fum CT	Aspergillus fumigatus Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5621475			Epicoccum	Epicoccum Purpurascens Allergen		
	15700-8	5621450	Epicoccum CL	Epicoccum Purpurascens Allergen	See below	See below
	6111-9	5621425	Epicoccum CT	Epicoccum Purpurascens Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5606175			Cladosporium	Cladosporium herbarum (Hormodendrum) Allergen		
	N/A	5606150	Cladosporium CL	Cladosporium herbarum (Hormodendrum) Allergen	See below	See below
	53760-5	5606125	Cladosporium CT	Cladosporium herbarum (Hormodendrum) Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5621325			Mucor Racemosus	Mucor Racemosus Allergen		
	15862-6	5621300	Mucor Race CL	Mucor Racemosus Allergen	See below	See below
	6182-0	5621275	Mucor Race CT	Mucor Racemosus Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
Remove component	5616830		Fusarium Moniliforme	Fusarium Moniliforme Allergen		
	15722-2	5616820	Fusarium CL	Fusarium Moniliforme Allergen	See below	See below
	6121-8	5616810	Fusarium CT	Fusarium Moniliforme Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
Remove component	5610775		Curvularia Lunata	Curvularia Lunata Allergen		
	15670-3	5610750	Curvularia CL	Curvularia Lunata Allergen	See below	See below
	42198-2	5610725	Curvularia CT	Curvularia Lunata Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
Remove component	5616780		Aur Pullulans	Aureobasidium Pullulans Allergen		
	15964-0	5616770	Aur Pullulan CL	Aureobasidium Pullulans Allergen	See below	See below
	6029-3	5616760	Aur Pullulan CT	Aureobasidium Pullulans Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
6906887			Cephalosporium	Cephalosporium		
	6066-5	6906887	Cephalosporium	Cephalosporium IgE	See below	See below
6906981			Corn Smut	Corn Smut IgE		
	N/A	6906983	Corn Smut CT	Corn Smut Count	See below	See below
	N/A	6906989	Corn Smut CL	Corn Smut Class	See below	See below
5622400			Pen Chrysogenum	Penicillium Chrysogenum Allergen		
	15924-4	5622375	Penicillium CL	Penicillium Chrysogenum Allergen	See below	See below
	6212-5	5622350	Penicillium CT	Penicillium Chrysogenum Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
Add component	6907249		Chaetomium	Chaetomium globosum IgE		
	7200-9	6907249	Chaetomium	Chaetomium IgE	See below	See below
Add component	2020850		IgE	Immunoglobulin IgE		
	N/A	2023100	IgE	IgE	See below	See below
	N/A	2023125	IgE Interp	IgE Interpretation	See below	See below
Add component	6907635		Bipol spicifera	Bipolaris spicifera IgE		
	42198-2		Bipol spicifera	Bipolaris spicifera IgE	See below	See below

Reference Ranges with Units of Measure are available upon request

Modify Test		Effective : Immediately		- Please update your EMR catalog with those appropriate to your practice -									
Test Code	Test Name	Mnemonic	Category/Type	Careset Component Information									
6907457	ENT Allergy Panel 2 (Moderate panel)	ENT PNL 2	Careset	Order Code	LOINC	Result Code	Mnemonic	Test Name	Reference Range	UOM	Result Type		
3 careset components added 4 careset components removed No other test changes				5606795			Elm Tree	Elm Tree Allergen					
					15697-6	5506785	Elm Tree CL	Elm Tree Allergen	See below		See below		
					6109-3	5606775	Elm Tree CT	Elm Tree Count	See below		See below		
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below		
					5613000		Silver Birch Tree	Common Silver Birch Tree Allergen					
					15283-5	5612075	Silver Birch CL	Common Silver Birch Tree Allergen	See below		See below		
					15579-6	5612050	Silver Birch CT	Common Silver Birch Tree Count	See below		See below		
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below		
					5611750		Box Elder	Box Elder (Maple) Tree Allergen					
					15585-3	5611725	Box Elder CL	Box Elder Tree Allergen	See below		See below		
					7155-5	5611700	Box Elder CT	Box Elder Tree Count	See below		See below		
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below		
					5607475		Cottonwood Tree	Cottonwood Tree Allergen					
					15659-6	5607450	Cottonwood CL	Cottonwood Tree Allergen	See below		See below		
					6090-5	5607425	Cottonwood CT	Cottonwood Tree Count	See below		See below		
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below		
					5518175		Mulberry Tree Red	Mulberry Tree Red Allergen					
					21400-7	5518150	Mulberry Red CL	Mulberry Tree Red Allergen	See below		See below		
					7514-3	5528650	Mulberry Red CT	Mulberry Tree Red Count	See below		See below		
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below		
					5528010		Mulberry Tree White	Mulberry Tree White Allergen					
					N/A	5528005	Mulberry Wht CT	Mulberry Tree White Allergen	See below		See below		
					N/A	5528015	Mulberry Wht CL	Mulberry Tree White Count	See below		See below		
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below		
					5611850		Mountain Juniper	Mountain Juniper Tree Allergen					
					15615-8	5611825	Mtn Juniper CL	Mountain Juniper Tree Allergen	See below		See below		
					6178-8	5611800	Mtn Juniper CT	Mountain Juniper Tree Count	See below		See below		
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below		
					5606750		Pecan, Hickory	Pecan, Hickory Tree Allergen					
					15774-3	5606725	Pecan Tr CL	Pecan, Hickory Tree Allergen	See below		See below		
					6209-1	5606700	Pecan Tr CT	Pecan, Hickory Tree Count	See below		See below		
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below		
					5621250		Privet Tree	Privet Tree Allergen					
					15960-8	5619050	Privet Tree CL	Privet Tree Allergen	See below		See below		
					7632-3	5619025	Privet Tree CT	Privet Tree Count	See below		See below		
					8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below		
					5613519		White Pine	White Pine Tree IgE					
					6282-8	5613519	White Pine Tr IgE	White Pine Tr IgE	See below		See below		
					5613225		Mesquite Tree	Mesquite Tree Allergen					
					15843-6	5613200	Mesquite Tr CL	Mesquite Tree Allergen	See below		See below		
		6173-9	5613175	Mesquite Tr CT	Mesquite Tree Count	See below		See below					
		8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below					
		5606325		Bermuda Grass	Bermuda Grass Allergen								
		15738-8	5606300	Bermuda CL	Bermuda Grass Allergen	See below		See below					
		6041-8	5606275	Bermuda CT	Bermuda Grass Count	See below		See below					
		8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below					
		5605965		Johnson Grass	Johnson Grass Allergen								
		15743-8	5605950	Johnson Gr CL	Johnson Grass Allergen	See below		See below					
		6152-3	5605925	Johnson Gr CT	Johnson Grass Count	See below		See below					
		8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below					
		5611950		Timothy Grass	Timothy Grass Allergen								
		16054-9	5611925	Timothy Gr CL	Timothy Grass Allergen	See below		See below					
		6265-3	5611900	Timothy Gr CT	Timothy Grass Count	See below		See below					
		8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below					
		5609225		Cocklebur	Cocklebur Allergen								
		15646-3	5609200	Cocklebur CL	Cocklebur Allergen	See below		See below					
		6077-2	5609175	Cocklebur CT	Cocklebur Count	See below		See below					
		8251-1	5606925	Allergen Interp	Allergen Interpretation	See below		See below					
		5582300		Dandelion	Dandelion Allergen								
		15676-0	5582250	Dandelion CL	Dandelion Allergen	See below		See below					
		6097-0	5582200	Dandelion CT	Dandelion Count	See below		See below					

6907457 ENT Allergy Panel 2 (Moderate panel) --- Continued on next page

6907457 ENT Allergy Panel 2 (Moderate panel) --- Continued from previous page

5607025			Giant Ragweed	Giant Ragweed Allergen		
	15978-0	5607000	Giant Ragweed CL	Giant Ragweed Allergen	See below	See below
	6124-2	5606975	Giant Ragweed CT	Giant Ragweed Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5621200			Firebush (Kochia)	Firebush (Kochia) Allergen		
	N/A	5621225	Firebush CL	Firebush (Kochia) Allergen	See below	See below
	7491-4	5621175	Firebush CT	Firebush (Kochia) Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5609300			Goosefoot Lamb's Quarte	Lamb's Quarters Goosefoot Allergen		
	15805-5	5609275	Lamb Quarter CL	Lamb's Quarters Goosefoot Allergen	See below	See below
	6156-4	5609250	Lamb Quarter CT	Lamb's Quarters Goosefoot Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5606835			Rough Marshelder	Rough Marshelder Allergen		
	15693-5	5606825	Rough Marsh CL	Rough Marshelder Allergen	See below	See below
	6232-3	5606800	Rough Marsh CT	Rough Marshelder Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5644025			Common Pigweed	Common Pigweed Allergen		
	15936-8	5644075	Pigweed CL	Common Pigweed Allergen	See below	See below
	7604-2	5644050	Pigweed CT	Common Pigweed Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5607625			Plantain (English)	Plantain (English), Ribwort Allergen		
	15952-5	5607600	Plantain Rib CL	Plantain (English), Ribwort Allergen	See below	See below
	6110-1	5607575	Plantain Rib CT	Plantain (English), Ribwort Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5607100			Russian Thistle	Saltwort, Russian Thistle Allergen		
	16052-3	5607072	Russian This CL	Saltwort, Russian Thistle Allergen	See below	See below
	6234-9	5607050	Russian This CT	Saltwort, Russian Thistle Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5604625			Common Ragweed	Common Ragweed Allergen		
	15975-6	5604600	Ragweed CL	Common Ragweed Allergen	See below	See below
	6085-5	5604575	Ragweed CT	Common Ragweed Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5612025			Mugwort	Mugwort Allergen		
	15863-4	5612000	Mugwort CL	Mugwort Allergen	See below	See below
	6183-8	5611975	Mugwort CT	Mugwort Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5606250			Cat Dander	Cat Dander Allergen		
	15609-1	5606225	Cat Dander CL	Cat Dander Allergen	See below	See below
	6833-8	5606200	Cat Dander CT	Cat Dander Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5613349			Cockroach American	Cockroach American IgE		
	30170-5	5613349	Cockroach Amer IgE	Cockroach American IgE	See below	See below
5606025			German Cockroach	German Cockroach Allergen		
	15647-1	5606000	Gr Roach CL	German Cockroach Allergen	See below	See below
	24139-8	5605975	Gr Roach CT	German Cockroach Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5605800			Dog Dander	Dog Dander Allergen		
	15685-1	5605750	Dog Dander CL	Dog Dander Allergen	See below	See below
	6098-8	5605775	Dog Dander CT	Dog Dander Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5604985			Dust Mite (D. farinae)	Dust Mite (D. farinae) Allergen		
	15680-2	5604975	Dust Mite CL	House Dust Mite Allergen	See below	See below
	6095-4	5604950	Dust Mite CT	House Dust Mite Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5607550			House Dust Mite	House Dust Mite Allergen		
	15682-8	5607525	House Dust D Petrony M	House Dust Mite D. Petrony Class	See below	See below
	6096-2	5607500	House Dust D Petrony M	House Dust Mite D. Petrony Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5606550			Alternaria Alternata	Alternaria Alternata Allergen		
	15530-9	5606525	Alt Ten CL	Alternaria Alternata Allergen	See below	See below
	6020-2	5606500	Alt Ten CT	Alternaria Alternata Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5606400			Aspergillus fumigatus	Aspergillus fumigatus Allergen		
	15549-9	5606375	Asp Fum CL	Aspergillus fumigatus Allergen	See below	See below
	6025-1	5606350	Asp Fum CT	Aspergillus fumigatus Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below
5621475			Epicoccum	Epicoccum Purpurascens Allergen		
	15700-8	5621450	Epicoccum CL	Epicoccum Purpurascens Allergen	See below	See below
	6111-9	5621425	Epicoccum CT	Epicoccum Purpurascens Count	See below	See below
	8251-1	5606925	Allergen Interp	Allergen Interpretation	See below	See below

6907457 ENT Allergy Panel 2 (Moderate panel) --- Continued on next page

6907457 ENT Allergy Panel 2 (Moderate panel) --- Continued from previous page

	5606175			Cladosporium	Cladosporium herbarum (Hormodendrum) Allergen		
	N/A	5606150		Cladosporium CL	Cladosporium herbarum (Hormodendrum) Allergen	See below	See below
	53760-5	5606125		Cladosporium CT	Cladosporium herbarum (Hormodendrum) Count	See below	See below
	8251-1	5606925		Allergen Interp	Allergen Interpretation	See below	See below
Remove component	5621325			Mucor Racemosus	Mucor Racemosus Allergen		
	15862-6	5621300		Mucor Race CL	Mucor Racemosus Allergen	See below	See below
	6182-0	5621275		Mucor Race CT	Mucor Racemosus Count	See below	See below
	8251-1	5606925		Allergen Interp	Allergen Interpretation	See below	See below
Remove component	5616830			Fusarium Moniliforme	Fusarium Moniliforme Allergen		
	15722-2	5616820		Fusarium CL	Fusarium Moniliforme Allergen	See below	See below
	6121-8	5616810		Fusarium CT	Fusarium Moniliforme Count	See below	See below
	8251-1	5606925		Allergen Interp	Allergen Interpretation	See below	See below
Remove component	5610775			Curvularia Lunata	Curvularia Lunata Allergen		
	15670-3	5610750		Curvularia CL	Curvularia Lunata Allergen	See below	See below
	42198-2	5610725		Curvularia CT	Curvularia Lunata Count	See below	See below
	8251-1	5606925		Allergen Interp	Allergen Interpretation	See below	See below
Remove component	5616780			Aur Pullulans	Aureobasidium Pullulans Allergen		
	15964-0	5616770		Aur Pullulan CL	Aureobasidium Pullulans Allergen	See below	See below
	6029-3	5616760		Aur Pullulan CT	Aureobasidium Pullulans Count	See below	See below
	8251-1	5606925		Allergen Interp	Allergen Interpretation	See below	See below
	6906887			Cephalosporium	Cephalosporium		
	6066-5	6906887		Cephalosporium	Cephalosporium IgE	See below	See below
	6906981			Corn Smut	Corn Smut IgE		
	N/A	6906983		Corn Smut CT	Corn Smut Count	See below	See below
	N/A	6906989		Corn Smut CL	Corn Smut Class	See below	See below
	5622400			Pen Chrysogenum	Penicillium Chrysogenum Allergen		
	15924-4	5622375		Penicillium CL	Penicillium Chrysogenum Allergen	See below	See below
	6212-5	5622350		Penicillium CT	Penicillium Chrysogenum Count	See below	See below
	8251-1	5606925		Allergen Interp	Allergen Interpretation	See below	See below
Add component	6907249			Chaetomium	Chaetomium globosum IgE		
	7200-9	6907249		Chaetomium	Chaetomium IgE	See below	See below
Add component	2020850			IgE	Immunoglobulin IgE		
	N/A	2023100		IgE	IgE	See below	See below
	N/A	2023125		IgE Interp	IgE Interpretation	See below	See below
Add component	6907635			Bipol spicifera	Bipolaris spicifera IgE		
	42198-2			Bipol spicifera	Bipolaris spicifera IgE	See below	See below

Reference Ranges with Units of Measure are available upon request

Modify Test Effective : Immediately - Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
3607775	Cholinesterase, Serum	CHOLINES S	Group						
Change in the Result Code from 3607775 to 3604200 No other test changes				2098-2	3607775	Cholinesterase	Cholinesterase	by report IU/L	Numeric XXXXXXXX
				2098-2	3604200	Pseudocholin	Pseudocholinesterase	by report IU/L	Numeric XXXXXXXX

Modify Test Effective : 11/16/2020 - Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	Careset Component Information						
Order Code	LOINC	Result Code	Mnemonic	Test Name	Reference Range UOM	Result Type				
6907567	Respiratory Virus Mini-Panel by PCR	Mini Vresp PCR		Careset Component Information						
One careset component added No other test changes				5565570			FLU A PCR	Influenza A by PCR		
					34487-9	5565570	FLU A PCR	Influenza A by PCR	See below	See below
				5565580			FLU B PCR	Influenza B by PCR		
					40982-1	5565580	FLU B PCR	Influenza B by PCR	See below	See below
				5565560			RSV PCR	RSV by PCR		
					40988-8	5565560	RSV PCR	RSV by PCR	See below	See below
				5575875			Metapneumo	Metapneumovirus PCR		
					38917-1	5575875	Metapneumo PCR	Metapneumovirus PCR	See below	See below
Add component	5565555			Adenov PCR			Adenov PCR	Adenovirus by PCR		
					39528-5	5565555	Adenov PCR	Adenovirus PCR	See below	See below

Reference Ranges with Units of Measure are available upon request

Modify Test Effective : 11/16/2020 - Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6906659	Adenosine Deaminase, Peritoneal Fluid	Adeno D Periton	Group						
Prompt has been discontinued No other test changes				49759-4	6906659	Adeno D Periton	Adenosine Deaminase, Peritoneal Fluid	X.X u/L	Free text
				N/A	6906689	Fluid-Type	Fluid-Type	N/A	Free text

Modify Test		Effective : 11/16/2020		- Please update your EMR catalog with those appropriate to your practice -						
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
3660750	Adenosine Deaminase, Pleural Fluid	ADENO DEAM	Group	35704-6	3660750	Adeno D Pleural	Adenosine Deaminase, Pleural Fluid	X.X u/L	Free text	
Prompt has been discontinued No other test changes				Discontinue	N/A	6906689	Fluid-Type	Fluid-Type	N/A	Free text

Modify Test		Effective : 11/16/2020		- Please update your EMR catalog with those appropriate to your practice -						
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
6904899	Adenosine Deaminase, Pericardial Fluid	Adeno D Pericard	Group	49760-2	49760-2	Adeno D Pericard	Adenosine Deaminase, Pericardial Fluid	X.X u/L	Free text	
Prompt has been discontinued No other test changes				Discontinue	N/A	6906689	Fluid-Type	Fluid-Type	N/A	Free text

Modify Test		Effective : 11/16/2020		- Please update your EMR catalog with those appropriate to your practice -						
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
5592675	Liver Cytosol (LC-1) Autoantibody	Liver Cyto	Detail	13175-5	5592675	Liver Cyto	Liver Cytosol Autoantibody	N/A	Alpha	
CPT has been updated No other test changes				CPT: 84182						

Modify Test		Effective : 11/16/2020		- Please update your EMR catalog with those appropriate to your practice -						
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
6906191	Myositis Extended Panel	Myositis Panel	Group	31625-7	5513216	SSA 52 (Ro) IgG	SSA 52 (Ro) ENA Ab IgG	See below	See below	
CPT: 83516 x 12 86235 x 6 CPT has been updated No other test changes				53019-6	5513215	SSA 60 (Ro) IgG	SSA 60 (Ro) Ab IgG	See below	See below	
				38369-5	5513217	Ribonucleic P Ab	Ribonucleic P Ab	See below	See below	
				14235-6	5513220	JO-1 Auto Ab	JO-1 Autoantibody	See below	See below	
				33771-7	5513140	PL-12 Auto Ab	PL-12 Autoantibody	See below	See below	
				33772-5	5513130	PL-7 Auto Ab	PL-7 Autoantibody	See below	See below	
				45152-6	5513190	OJ Auto Ab	OJ Autoantibody	See below	See below	
				45149-2	5513180	EJ Auto Ab	EJ Autoantibody	See below	See below	
				33921-8	5513210	SRP Auto Ab	Signal Recognition Particle (SRP) Autoantibody	See below	See below	
				18484-6	5513170	Ku Auto Ab	Ku Autoantibody	See below	See below	
				81732-0	3806050	Pm-1 Ab	PM-SCL Antibody	See below	See below	
				49963-2	6906193	Fibrillarin IgG	Fibrillarin (U3 RNP) Ab, IgG	See below	See below	
				18485-3	5513160	Mi-2 Auto Ab	Mi-2 Autoantibody	See below	See below	
				N/A	5513211	P155/140 Ab	P155/140 Ab	See below	See below	
				N/A	6906195	TIF1-gamma Ab	TIF1-gamma Antibody	See below	See below	
				N/A	6906197	SAE1 Antibody	SAE1 (SUMO activating enzyme) Antibody	See below	See below	
				N/A	6906199	MDA5 Antibody	MDA5 (CADM-140) Antibody	See below	See below	
N/A	6906201	NXP-2 Ab	NXP-2 (Nuclear matrix protein-2) Ab	See below	See below					
48767-8	6906203	Myositis Pn Int	Myositis Panel Interpretive Data	See below	See below					
Reference Ranges with Units of Measure are available upon request										

Modify Test		Effective : 11/16/2020		- Please update your EMR catalog with those appropriate to your practice -						
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
6904901	Bupropion (Wellbutrin, Zyban) S/P	Bupropion S/P	Group	6706-6	6904911	Bupropion Lvl	Bupropion Level	XXXX ng/mL	Free text	
One new result component No other test changes				Add component	9418-5	6907313	Hydroxybup Lvl	Hydroxybupropion Level	XXXX ng/mL	Free text

INACTIVATED TESTS - Please update your EMR catalog with those appropriate to your practice

INACTIVE		Effective : 11/16/2020		Replacement Information					
Inactive	Inactive	Inactive		Test Code	Mnemonic	Test Name			
4500450	Cortisol AM	CORTISL A		4400350	TOT CORTIS	Cortisol Serum Total			

The Cortisol AM and Cortisol PM tests will be discontinued and replaced by the Cortisol Total assay.
The Cortisol Total will be reported with a footnote containing the reference ranges for AM and PM collection.

INACTIVE Effective : 11/16/2020						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
4500500	Cortisol PM	CORTISL P	4400350	TOT CORTIS	Cortisol Serum Total	
<p>The Cortisol AM and Cortisol PM tests will be discontinued and replaced by the Cortisol Total assay. The Cortisol Total will be reported with a footnote containing the reference ranges for AM and PM collection.</p>						

INACTIVE Effective : 11/16/2020						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6906789	Breast & Ovarian Hearedity Cancer Syndrome (BRCA1/BRCA2) Sequencing	BRCA1/2 Seq	6907511	BRCA1/2 NGS S/D/D	BRCA1/2 NGS Sequencing and Del/Dup	
<p>See New Test section for test build information</p>						

INACTIVE Effective : 11/16/2020						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6906815	Breast & Ovarian Hearedity Cancer Syndrome (BRCA1/BRCA2) Sequencing & Deletion/Duplication	BRCA1/2 Seq/Del	6907511	BRCA1/2 NGS S/D/D	BRCA1/2 NGS Sequencing and Del/Dup	
<p>See New Test section for test build information</p>						

INACTIVE Effective : 11/16/2020						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6004175	Coccidioides Antibody, Immunodiffusion	Coccid Ab ID	6907499	Coccid Ab ID	Coccidioides Ab by ID	
<p>See New Test section for test build information</p>						