

NEW TESTS - Please update your EMR catalog with those appropriate to your practice
New Test Effective : 4/18/2022

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6907558	COV2, Flu A/B, RSV Respiratory 4Plex	COV2 FLU RSV 4P	Group						
CPT: 87637									
1 Specimen sites: Nasopharyngeal (NP), Oropharyngeal (OP), Deep Nasal/Mid-turbinate (MT) or Anterior Nares (AN) swab. If collecting MT or AN specimens, both nares should be generously sampled with a single swab and placed in a single tube to maximize virus particle collection.									
2 Any swab of appropriate size and configuration with a synthetic or rayon tip and plastic shaft can be used.									
3 Do NOT submit swabs with calcium alginate tips or swabs with wood shafts. All Other Swab types: Suitable for all specimen sites listed, depending on size of swab Refrigerate specimens if they will be received in lab within 72 hours of collection. If there will be delayed transport, freeze specimen at -20°C									
4 Media types: VTM, UTM, M6, or sterile saline (ESwab collection kits should no longer be used for specimen collection for the COV2, Flu A/B, RSV Respiratory 4Plex assay. The small volume of this collection device can cause technical issues with PCR and precludes any repeat or additional testing. Contact your RML representative for alternate collection devices.)									
5 Swabs must be placed in 1-3 mL of liquid media, broken or cut at the top of the tube, and tubes tightly capped. Specimen tubes should be labeled with two unique patient identifiers. requisition or place electronic order. Place the tube containing the swab in a biohazard bag and submit with paper									
6 Please keep specimen refrigerated until pick-up.									
7 Causes for rejection: Swabs with calcium alginate tips; swabs with wooden shafts; refrigerated samples greater than 72 hours old; room temperature swabs greater than 24 hours old; improperly labeled; grossly contaminated; broken or leaking tube; collection with substances inhibitory to PCR including heparin, hemoglobin, ethanol, EDTA concentrations >0.01M; unacceptable transport media.									
8 Couriers will transport to Regional Medical Core Laboratory. Call RML at 918-744-2232 to request specimen pickup if not routinely scheduled at your location.									

New Test Effective : 5/16/2022

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6907547	Immune Complex, C3 Fragments	C3 Fragments	Detail	10864-7	6907547	C3 Fragments	Immune Complex, C3 Fragments	<15 ugEq/mL	Numeric XXXXX

CPT: 86332 Collect: Plain red or Serum separator tube (SST)
 Specimen Prep: Allow complete clotting of red blood cells (up to 1 hour), then separate serum from cells within 30 minutes and freeze immediately.
 Transfer 1mL (0.5mL) of serum to a pour off tube.
 If ordered in conjunction with a C1q Binding Assay, submit 2 separate 1mL aliquots of serum.
CRITICAL FROZEN! Separate specimens must be submitted when multiple tests are ordered.
 Unacceptable: Non-frozen specimens. Specimens exposed to repeated freeze/thaw cycles. Grossly hemolyzed, lipemic and icteric specimens.
 Stability: Ambient-Unacceptable; Refrigerated-Unacceptable; Frozen- 20 days
 Method: Semi-Quant ELISA
 Set-up: Sat, Mon, Thur TAT: 2-9 days after set-up

New Test Effective : 5/16/2022									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6908321	Drug Meconium Qual Non-COC	Drug M Qual NCOC	Group						
CPT: 80326, 80347, 80364, 80355 (Alt:G0481) Collect: All meconium (blackish material) excreted until milk/formula-based stool (yellow-green) appears. Specimen Prep: Transport all available meconium (4g is preferred) to routine urine collection cup or Security Kit for Meconium/Umbilical Drug Detection. Transport Temp: Refrigerated Unacceptable: Unknown fluids, pharmaceutical preparation, and breast milk. Diapers, cotton swabs, baby wipes, tongue depressors, bottles. Stability: Ambient:1 week; Refrigerated: 3 weeks; Frozen: 1 year Method: Qualitative LC/MS/MS Set-up: Sun-Sat TAT: 1-3 days after set-up				32099-4	6908323	Mec 6-Acetylmor	6-Acetylmorphine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				61031-1	6908325	Mec 7A-clonazep	7-Aminoclonazepam, MEC, Qual	< 5 ng/g	Numeric XXXXXXXX
				61037-8	6908327	Mec AOH-Alpraz	Alpha-OH-Alprazolam, MEC, Qual	< 5 ng/g	Numeric XXXXXXXX
				N/A	6908329	Mec AOH-Midazo	Alpha-OH-Midazolam, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				61038-6	6908331	Mec Alprazolam	Alprazolam, MEC, Qual	< 5 ng/g	Numeric XXXXXXXX
				8144-8	6908333	Mec Amphetamine	Amphetamine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				8187-7	6908335	Mec Benzoylcego	Benzoylcegonine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				N/A	6908337	Mec Buprenorph	Buprenorphine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				32057-2	6908339	Mec Butalbital	Butalbital, MEC, Qual	< 50 ng/g	Numeric XXXXXXXX
				61039-4	6908341	Mec Clonazepam	Clonazepam, MEC, Qual	< 5 ng/g	Numeric XXXXXXXX
				57945-8	6908343	Mec Cocaethylen	Cocaethylene, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				40527-4	6906305	Mec Cocaine	Cocaine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				29278-9	6908347	Mec Codeine	Codeine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				61074-1	6908349	Mec Diazepam	Diazepam, MEC, Qual	< 5 ng/g	Numeric XXXXXXXX
				N/A	6908351	Mec Dihydrocod	Dihydrocodeine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				N/A	6908353	Mec Methad Met	Methadone Metabolite, MEC, Qual	< 10 ng/g	Numeric XXXXXXXX
				61042-8	6908355	Mec Fentanyl	Fentanyl, MEC, Qual	< 10 ng/g	Numeric XXXXXXXX
				N/A	6908357	Mec Gabapentin	Gabapentin, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				57981-3	6908359	Mec Hydrocodone	Hydrocodone, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				57982-1	6908361	Mec Hydromorph	Hydromorphone, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				61044-4	6908363	Mec Lorazepam	Lorazepam, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				40481-4	6908365	Mec MDMA	MDMA- Ecstasy, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				75647-8	6908367	Mec Meperidine	Meperidine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				43595-8	6906313	Mec Methadone	Methadone, MEC, Qual	< 10 ng/g	Numeric XXXXXXXX
				29286-2	6908371	Mec Methamphet	Methamphetamine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				61050-1	6908373	Mec Methylphen	Methylphenidate, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				N/A	6908375	Mec Midazolam	Midazolam, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				N/A	6908377	Mec m-OH-Benzoy	m-OH-Benzoylcegonine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				29287-0	6908379	Mec Morphine	Morphine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				N/A	6908381	Mec N-Tramadol	N-desmethyltramadol, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				N/A	6908383	Mec Naloxone	Naloxone, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				N/A	6908385	Mec Norbuprenor	Norbuprenorphine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				61051-9	6908387	Mec Nordiazepam	Nordiazepam, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				N/A	6908389	Mec Norhydrocod	Norhydrocodone, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				N/A	6908391	Mec Noroxycodon	Noroxycodone, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				N/A	6908393	Mec O-Tramadol	O-desmethyltramadol, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				61055-0	6908395	Mec Oxazepam	Oxazepam, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				32101-8	6908397	Mec Oxycodone	Oxycodone, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				N/A	6908399	Mec Oxymorphone	Oxymorphone, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX
				8232-1	6906309	Mec Phencyclid	Phencyclidine- PCP, MEC, Qual	< 10 ng/g	Numeric XXXXXXXX
			32108-3	6908403	Mec Phenobarbit	Phenobarbital, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX	
			N/A	6908405	Mec Phentermine	Phentermine, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX	
			N/A	6908407	Mec Tapentadol	Tapentadol, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX	
			61061-8	6908409	Mec Temazepam	Temazepam, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX	
			N/A	6908411	Mec Tramadol	Tramadol, MEC, Qual	< 20 ng/g	Numeric XXXXXXXX	
			N/A	6908413	Mec Zolpidem	Zolpidem, MEC, Qual	< 10 ng/g	Numeric XXXXXXXX	
			11502-2	6908415	EER Mec Report	EER Drug Detection Pan, MEC	By Report	Free text	
			58738-6	6908417	Drug Det Pnl Mec	Drug Detection Panel, MEC	By Report	Free text	

MODIFICATION TESTS - Please update your EMR catalog with those appropriate to your practice
Modify Test Effective : Immediately

- Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
5565570	Influenza A by PCR	FLU A PCR	Detail	34487-9	5565570	Flu A PCR	Influenza A by PCR	Not Detected	Alpha
Changing CPT code from 87798 to 87501 No other changes									

Modify Test Effective : Immediately

- Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
5565580	Influenza B by PCR	FLU B PCR	Detail	40982-1	5565580	Flu B PCR	Influenza B by PCR	Not Detected	Alpha
Changing CPT code from 87798 to 87501 No other changes									

Modify Test Effective : Immediately

- Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6906339	Lead Venous	Lead Venous	Detail	5671-3	3601650	Lead Lvl	Lead Level	<3.5 mcg/dL	Numeric
Reference range has been updated No other changes									

Modify Test Effective : Immediately

- Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6906337	Lead Capillary	Lead Capillary	Detail	10368-9	6906337	Lead Capillary	Lead Capillary	<3.5 mcg/dL	Numeric
Reference range has been updated No other changes									

Modify Test Effective : 5/16/2022

- Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
3600160	Angiotensin Converting Enzyme	ANGIOTEN	Detail	2742-5	3600160	ACE	Angiotensin Converting Enzyme	See Below	Numeric
Reference range has been updated No other changes									
									>= 18 yrs: 16-85 U/L

Modify Test Effective : 05/16/2022 - Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type																																																								
5572555	Collagen Type I C-Telopeptide(CTx)	C TEOPEP	Detail	41171-0	5572555	C Telo pep	C Telopeptide b-crosslinked	See Below	Numeric XXXXXXXX																																																								
<p>Reference range has been updated No other changes</p> <table border="0"> <tr> <td></td> <td>Age</td> <td>Female</td> <td>Male</td> </tr> <tr> <td></td> <td>6 months-6 years</td> <td>500-1800 pg/mL</td> <td>238-1019 pg/mL</td> </tr> <tr> <td></td> <td>7-9 years</td> <td>566-1690 pg/mL</td> <td>238-1019 pg/mL</td> </tr> <tr> <td></td> <td>10-12 years</td> <td>503-2077 pg/mL</td> <td>238-1019 pg/mL</td> </tr> <tr> <td></td> <td>13-15 years</td> <td>160-1590 pg/mL</td> <td>238-1019 pg/mL</td> </tr> <tr> <td></td> <td>16-17 years</td> <td>167-933 pg/mL</td> <td>238-1019 pg/mL</td> </tr> <tr> <td></td> <td>18-29 years</td> <td></td> <td>238-1019 pg/mL</td> </tr> <tr> <td></td> <td>30-39 years</td> <td></td> <td>225-936 pg/mL</td> </tr> <tr> <td></td> <td>40-49 years</td> <td></td> <td>182-801 pg/mL</td> </tr> <tr> <td></td> <td>50-59 years</td> <td></td> <td>161-737 pg/mL</td> </tr> <tr> <td></td> <td>60-69 years</td> <td></td> <td>132-752 pg/mL</td> </tr> <tr> <td></td> <td>>= 70 years</td> <td></td> <td>118-776 pg/mL</td> </tr> <tr> <td></td> <td>Premenopausal</td> <td>136-689 pg/mL</td> <td></td> </tr> <tr> <td></td> <td>Postmenopausal</td> <td>177-1015 pg/mL</td> <td></td> </tr> </table>											Age	Female	Male		6 months-6 years	500-1800 pg/mL	238-1019 pg/mL		7-9 years	566-1690 pg/mL	238-1019 pg/mL		10-12 years	503-2077 pg/mL	238-1019 pg/mL		13-15 years	160-1590 pg/mL	238-1019 pg/mL		16-17 years	167-933 pg/mL	238-1019 pg/mL		18-29 years		238-1019 pg/mL		30-39 years		225-936 pg/mL		40-49 years		182-801 pg/mL		50-59 years		161-737 pg/mL		60-69 years		132-752 pg/mL		>= 70 years		118-776 pg/mL		Premenopausal	136-689 pg/mL			Postmenopausal	177-1015 pg/mL	
	Age	Female	Male																																																														
	6 months-6 years	500-1800 pg/mL	238-1019 pg/mL																																																														
	7-9 years	566-1690 pg/mL	238-1019 pg/mL																																																														
	10-12 years	503-2077 pg/mL	238-1019 pg/mL																																																														
	13-15 years	160-1590 pg/mL	238-1019 pg/mL																																																														
	16-17 years	167-933 pg/mL	238-1019 pg/mL																																																														
	18-29 years		238-1019 pg/mL																																																														
	30-39 years		225-936 pg/mL																																																														
	40-49 years		182-801 pg/mL																																																														
	50-59 years		161-737 pg/mL																																																														
	60-69 years		132-752 pg/mL																																																														
	>= 70 years		118-776 pg/mL																																																														
	Premenopausal	136-689 pg/mL																																																															
	Postmenopausal	177-1015 pg/mL																																																															

Modify Test Effective : 05/16/2022 - Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type																					
6905139	S-100B Protein, Serum	S-100B Protein, Serum	Detail	47275-3	6905139	S100B Prtn Srm	S-100B Protein, Serum	See Below	Numeric XXXXXXXX																					
<p>Reference range has been updated No other changes</p> <table border="0"> <tr> <td></td> <td>Age</td> <td>Values</td> </tr> <tr> <td></td> <td>0 days-5 months</td> <td><=510 ng/mL</td> </tr> <tr> <td></td> <td>5-10 months</td> <td><=3350 ng/mL</td> </tr> <tr> <td></td> <td>10-25 months</td> <td><=230 ng/mL</td> </tr> <tr> <td></td> <td>25 months-3 years</td> <td><=170 ng/L</td> </tr> <tr> <td></td> <td>3-19 years</td> <td><=160 ng/L</td> </tr> <tr> <td></td> <td>>=19 years</td> <td>0-96 ng/L</td> </tr> </table>											Age	Values		0 days-5 months	<=510 ng/mL		5-10 months	<=3350 ng/mL		10-25 months	<=230 ng/mL		25 months-3 years	<=170 ng/L		3-19 years	<=160 ng/L		>=19 years	0-96 ng/L
	Age	Values																												
	0 days-5 months	<=510 ng/mL																												
	5-10 months	<=3350 ng/mL																												
	10-25 months	<=230 ng/mL																												
	25 months-3 years	<=170 ng/L																												
	3-19 years	<=160 ng/L																												
	>=19 years	0-96 ng/L																												

Modify Test Effective : 05/16/2022 - Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6906341	Umbilical Cord Tissue, Drug Screen Non-Chain of Custody, Qual	UCT Drug NCOC	Group						
<p>CPT codes have been updated: From 80307 to 80326, 80347, 80364, 80355 (Alt: G0481) No other changes</p> <p style="text-align: center;">For Result Code information contact RML Help Desk</p>									

Modify Test Effective : 05/16/2022 - Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
3488394	MEN2, RET Gene Sequencing	MEN2 Gn Seq	Group						
<p>CPT code has been updated: from 81405 to 81406 No other changes</p> <p style="text-align: center;">For Result Code information contact RML Help Desk</p>									

INACTIVATED TESTS - Please update your EMR catalog with those appropriate to your practice
INACTIVE Effective : Immediately

Inactive	Inactive	Inactive	Replacement Information		
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name
6905105	Periprosthetic Joint Infection (PJI) Detection (Synovasure or Alpha Defensin	SYNOVA PJI	N/A	N/A	N/A

There is no available replacement at this time

INACTIVE Effective : 5/16/2022

Inactive	Inactive	Inactive	Replacement Information		
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name
6905673	Connexin 26 (GJB2) Sequencing	Connexin 26 Seq	N/A	N/A	N/A

There is no available replacement at this time

INACTIVE Effective : 5/16/2022

Inactive	Inactive	Inactive	Replacement Information		
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name
6905097	Pancreatitis (SPINK1) Sequencing	SPINK1 FGS	N/A	N/A	N/A

There is no available replacement at this time

INACTIVE Effective : 5/16/2022

Inactive	Inactive	Inactive	Replacement Information		
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name
5500430	Raji Cell Immune Complex Assay	RAJI CELL	6907547	C3 Fragments	Immune Complex, C3 Fragments

See New Test section for test build information

INACTIVE Effective : 5/16/2022

Inactive	Inactive	Inactive	Replacement Information		
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name
6906301	Drug Screen Meconium, Non-COC	DRUG M NCOC	6908321	Drug M Qual NCOC	Drug Meconium Qual Non-COC

See New Test section for test build information

INACTIVE Effective : 5/16/2022

Inactive	Inactive	Inactive	Replacement Information		
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name
6906451	.Rflx Amphetamines Conf, Mec Bill	.RfCB Amphetamin	6908321	Drug M Qual NCOC	Drug Meconium Qual Non-COC

See New Test section for test build information

INACTIVE Effective : 5/16/2022						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6906453	.Rflx Barbiturates Conf, Mec Bill	.RfCB Barbiturat	6908321	Drug M Qual NCOC	Drug Meconium Qual Non-COC	
See New Test section for test build information						
INACTIVE Effective : 5/16/2022						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6906457	.Rflx Benzodiazepine Conf, Mec Bill	.RfCB Benzodiaze	6908321	Drug M Qual NCOC	Drug Meconium Qual Non-COC	
See New Test section for test build information						
INACTIVE Effective : 5/16/2022						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6906449	.Rflx Phencyclidine Conf, Mec Bill	.RfCB Phencyclid	6908321	Drug M Qual NCOC	Drug Meconium Qual Non-COC	
See New Test section for test build information						
INACTIVE Effective : 5/16/2022						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6906461	.Rflx Marijuana Conf Mec, Bill	.RfCB Marijuana	6908321	Drug M Qual NCOC	Drug Meconium Qual Non-COC	
See New Test section for test build information						
INACTIVE Effective : 5/16/2022						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6906445	.Rflx Cocaine Conf, Mec Bill	.RfCB Cocaine	6908321	Drug M Qual NCOC	Drug Meconium Qual Non-COC	
See New Test section for test build information						
INACTIVE Effective : 5/16/2022						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6906447	.Rflx Opiates Conf, Mec Bill	.RfCB Opiates	6908321	Drug M Qual NCOC	Drug Meconium Qual Non-COC	
See New Test section for test build information						
INACTIVE Effective : 5/16/2022						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6906455	.Rflx Methadone Conf, Mec Bill	.RfCB Methadone	6908321	Drug M Qual NCOC	Drug Meconium Qual Non-COC	
See New Test section for test build information						
INACTIVE Effective : 5/16/2022						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6906459	.Rflx Buprenorphine Confm, Mec - Bill	.RfCB Buprenorph	6908321	Drug M Qual NCOC	Drug Meconium Qual Non-COC	
See New Test section for test build information						