

MODIFICATION TESTS - Please update your EMR catalog with those appropriate to your practice

Modify Test Effective : Immediately - Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
4008600	Hypoglycemic Panel	HYPOGLYC P	Group	48329-7	4008850	Chlorpromazine	Chlorpropamide Level	No change	No change	
CPT: G0480 Collect: Plain red or gray (sodium fluoride/potassium oxalate) Specimen Prep: Separate from cells ASAP or within 2 hours of collection. Transfer 1mL (0.3mL) serum or plasma to a transport tube. Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered. Storage/transport: Frozen or Refrigerated Unacceptable: Separator tubes Stability: Ambient: 48 hours; Refrigerated: 28 days; Frozen: 24 months Method: HPLC-LC/MS/MS Set-up: Varies TAT: 4-7 days after set-up at testing lab.				21566-5	4008900	Tolazamide Lvl	Tolazamide Level	No change	No change	
				48327-1	4009000	Glyburide Lvl	Glyburide Level	No change	No change	
				Remove this result	NA	4009050	Acetohexam Lvl	Acetohexamide Level	No change	No change
				21567-3	4008950	Tolbutami Lvl	Tolbutamide Level	No change	No change	
				48326-3	4004100	Glipizide Lvl	Glipizide Level	No change	No change	
				48325-5	4004750	Glimepiride Lvl	Glimepiride Level	No change	No change	
				49487-2	4005450	Nateglinide Lvl	Nateglinide Level	No change	No change	
				48328-9	4005500	Repaglinide Lvl	Repaglinide Level	No change	No change	
				Add this result	NA	4005501	Rosiglitazn Lvl	Rosiglitazone Level	By Report	Alpha
				Add this result	NA	4005502	Pioglitazon Lvl	Pioglitazone Level	By Report	Alpha