

NEW TESTS - Please update your EMR catalog with those appropriate to your practice

New Test Effective : 6/5/2017				Corrected Result Codes for the Gelatin Allergen (Bovine)					
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
5613375	Gelatin Allergen (Bovine)	GELATIN IC	Group	n/a	5613350	Gelatin CL	Gelatin Allergen	<0.35 kU/L	Numeric XXX.XX
CPT: 86003	1 mL (0.1 mL) Serum Clot Activator SST (Red/Gray or Tiger Top) Refrigerated			n/a	5613325	Gelatin CT	Gelatin Count	Class 0	Alpha
Patient: \$30.05	Specimen Stability: refrigerated 4 weeks, frozen 4 weeks. Unacceptable: Room Temperature			n/a	5606925	Allergen Interp	Allergen Interpretation	By Report	Alpha
<p><i>This is the suggested replacement to the discontinued test below [5508725] This is Not a New Test but you may not have this built in your system</i></p>				<p><i>See the Corrected Result codes for the Gelatin Allergy test</i></p>					

New Test Effective : 06/12/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6006452	Stratify JCV Antibody w/Index and Reflex to Inhibition Assay	JCV Ab/Index	Group	70173-0	6006453	JCV Antibody	JCV Antibody	Negative	Alpha
CPT: 86711	Collect: Serum from Separator Tube (SST) or Plasma from Lavender (EDTA)				6006454	JCV Index	JCV Index	<0.20	Numeric X.XX
Patient: \$0.00	Separate serum or plasma from cells. Transport 0.5 mL(0.25 mL) Serum or Plasma at Room Temperature			If Reflexed		Reflex Order Code:	6006473	JCV Antibody Inhibition Assay	
If Reflexed	Unacceptable Conditions: Heparinized specimens.			LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
CPT: 86711	Stability: Ambient: 7 Days; Refrigerated: 14 Days; Frozen: 90 Days			N/A	6006473	JCV Inhb Ab	JCV Antibody Inhibition Assay	By report	By report
Patient: \$0.00	Pricing is currently \$0.00 per the performing laboratory								

MODIFICATION TESTS - Please update your EMR catalog with those appropriate to your practice

Modify Test Effective : 6/5/2017																					
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type												
5960165	Chlamydia trachomatis DNA Rectal/Throat/Eye (Aptima Probe)	REF CHLAM APTI	Detail	43304-5	5960165	Chlam Aptima	Chlamydia Trachomatis by Amplified Detection	Negative	Alpha												
CPT: 87491	<i>Name Changed : Use This order specifically those collection types.</i>				6904843	APTIMA Media Type	APTIMA Transport Media Type	Free Text													
Patient \$129.75	<i>This test is to be used only for Rectal/Throat/Eye specimen types in an APTIMA Unisex Swab Specimen Collection kit.</i>				6904845	Source A2	Specimen Source	Free Text													
<p>Unacceptable Conditions: Large white swab included in APTIMA Unisex Swab Specimen Collection kit is for preparatory cleaning of the endocervix and is unacceptable for testing. Specimens in any transport media other than indicated above. Specimens in swab transport media without a swab. Specimen source is required. Storage/Transport Temperature: Refrigerated. Stability: Swab in APTIMA Swab Specimen Transport Tube: Ambient: 2 months; Refrigerated: 2 months; Frozen: 1 year</p>				<p><i>For Urine, Vaginal Swab, Male Urethral Swab, or Cervical Brush in ThinPrep Pap See the appropriate test below.</i></p> <table border="1"> <thead> <tr> <th>Test Code</th> <th>Mnemonic</th> <th>Test Name</th> </tr> </thead> <tbody> <tr> <td>6910111</td> <td>CT Urine</td> <td>C.trachomatis Urine</td> </tr> <tr> <td>6910109</td> <td>CT Swab</td> <td>C.trachomatis Swab</td> </tr> <tr> <td>6910107</td> <td>CT PAP</td> <td>C.trachomatis PAP</td> </tr> </tbody> </table> <p><i>Urine Specimen Swab Specimen ThinPrep Specimen</i></p> <p><i>For Test Build Information See the 05/22/2017 LiveWire.</i></p>						Test Code	Mnemonic	Test Name	6910111	CT Urine	C.trachomatis Urine	6910109	CT Swab	C.trachomatis Swab	6910107	CT PAP	C.trachomatis PAP
Test Code	Mnemonic	Test Name																			
6910111	CT Urine	C.trachomatis Urine																			
6910109	CT Swab	C.trachomatis Swab																			
6910107	CT PAP	C.trachomatis PAP																			

Modify Test Effective : 6/5/2017																					
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type												
5960185	Neisseria Gonorrhoeae DNA Rectal/Throat/Eye (Aptima Probe)	REF GC APTIMA	Detail	24111-7	5960185	GC Aptima	Neisseria Gonorrhoeae by Amplified Detection	Negative	Alpha												
CPT: 87591	<i>Name Changed : Use This order specifically those collection types.</i>				6904843	APTIMA Media Type	APTIMA Transport Media Type	Free Text													
Patient \$129.75	<i>This test is to be used only for Rectal/Throat/Eye specimen types in an APTIMA Unisex Swab Specimen Collection kit.</i>				6904845	Source A2	Specimen Source	Free Text													
<p>Unacceptable Conditions: Large white swab included in APTIMA Unisex Swab Specimen Collection kit is for preparatory cleaning of the endocervix and is unacceptable for testing. Specimens in any transport media other than indicated above. Specimens in swab transport media without a swab. Specimen source is required. Storage/Transport Temperature: Refrigerated. Stability: Swab in APTIMA Swab Specimen Transport Tube: Ambient: 2 months; Refrigerated: 2 months; Frozen: 1 year</p>				<p><i>For Urine, Vaginal Swab, Male Urethral Swab, or Cervical Brush in ThinPrep Pap See the appropriate test below.</i></p> <table border="1"> <thead> <tr> <th>Test Code</th> <th>Mnemonic</th> <th>Test Name</th> </tr> </thead> <tbody> <tr> <td>6910117</td> <td>GC Urine</td> <td>N.gonorrhoeae Urine</td> </tr> <tr> <td>6910115</td> <td>GC Swab</td> <td>N.gonorrhoeae Swab</td> </tr> <tr> <td>6910113</td> <td>GC PAP</td> <td>N.gonorrhoeae PAP</td> </tr> </tbody> </table> <p><i>Urine Specimen Swab Specimen ThinPrep Specimen</i></p> <p><i>For Test Build Information See the 05/22/2017 LiveWire.</i></p>						Test Code	Mnemonic	Test Name	6910117	GC Urine	N.gonorrhoeae Urine	6910115	GC Swab	N.gonorrhoeae Swab	6910113	GC PAP	N.gonorrhoeae PAP
Test Code	Mnemonic	Test Name																			
6910117	GC Urine	N.gonorrhoeae Urine																			
6910115	GC Swab	N.gonorrhoeae Swab																			
6910113	GC PAP	N.gonorrhoeae PAP																			

Modify Test Effective : 6/5/2017																					
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type												
5560335	GC and Chlamydia Rectal/Throat/Eye (Aptima Probe)	REF GC/CHL APT	Detail	43304-5	5960165	Chlam Aptima	Chlamydia Trachomatis by Amplified Detection	Negative	Alpha												
CPT: 87491; 87591 Patient \$259.50				24111-7	5960185	GC Aptima	Neisseria Gonorrhoeae by Amplified Detection	Negative	Alpha												
<p><i>Name Changed : Use This order specifically those collection types.</i></p> <p><i>This test is to be used only for Rectal/Throat/Eye specimen types in an APTIMA Unisex Swab Specimen Collection kit.</i></p> <p>Unacceptable Conditions: Large white swab included in APTIMA Unisex Swab Specimen Collection kit is for preparatory cleaning of the endocervix and is unacceptable for testing. Specimens in any transport media other than indicated above. Specimens in swab transport media without a swab. Specimen source is required. Storage/Transport Temperature: Refrigerated. Stability: Swab in APTIMA Swab Specimen Transport Tube: Ambient: 2 months; Refrigerated: 2 months; Frozen: 1 year</p>				Prompt	6904843	APTIMA Media Type	APTIMA Transport Media Type	Free Text													
				Prompt	6904845	Source A2	Specimen Source	Free Text													
				<p><i>For Urine, Vaginal Swab, Male Urethral Swab, or Cervical Brush in ThinPrep Pap See the appropriate test below.</i></p> <table border="1"> <thead> <tr> <th>Test Code</th> <th>Mnemonic</th> <th>Test Name</th> </tr> </thead> <tbody> <tr> <td>6910129</td> <td>CT/GC Urine</td> <td>C.trachomatis,N.gonorrhoeae Urine</td> </tr> <tr> <td>6910127</td> <td>CT/GC Swab</td> <td>C.trachomatis,N.gonorrhoeae Swab</td> </tr> <tr> <td>6910125</td> <td>CT/GC PAP</td> <td>C.trachomatis,N.gonorrhoeae PAP</td> </tr> </tbody> </table> <p><i>For Test Build Information See the 05/22/2017 LiveWire.</i></p>						Test Code	Mnemonic	Test Name	6910129	CT/GC Urine	C.trachomatis,N.gonorrhoeae Urine	6910127	CT/GC Swab	C.trachomatis,N.gonorrhoeae Swab	6910125	CT/GC PAP	C.trachomatis,N.gonorrhoeae PAP
Test Code	Mnemonic	Test Name																			
6910129	CT/GC Urine	C.trachomatis,N.gonorrhoeae Urine																			
6910127	CT/GC Swab	C.trachomatis,N.gonorrhoeae Swab																			
6910125	CT/GC PAP	C.trachomatis,N.gonorrhoeae PAP																			
Modify Test Effective : 6/5/2017																					
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type												
6900571	Thiopurine Metabolites (6-TG, 6-MMP)	Thiopurine Met	Group	32660-3	6900573	6-TG Nucleotid	6-TG Nucleotides RBC	230-400 pmol/8 x 108 RBC	Numeric XXXXX												
CPT: 80299 Patient \$ 345.00				32654-6	6900575	6-MMP Nucleotid	6-MMP Nucleotides RBC	Less than 5701 pmol/8 x 108 RBC	Numeric XXXXX												
<p><i>Change in Specimen Requirements and TAT</i></p> <p><i>Collection: 5.0 mL (2.5mL) Whole Blood. Refrigerated</i></p> <p><i>Specimen: EDTA (Lavender Top)</i></p> <p><i>Stability: Room temperature: 24 hours, Refrigerated: 5 days, Frozen: Unacceptable</i></p> <p><i>Set-up: Varies TAT: 3-8 days</i></p>																					
Modify Test Effective : 6/5/2017																					
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type												
5577425	Platelet Refractory Antibody	PLT REFAB	Detail	5577425		PLT Ref Ab	Platelet Refractory Antibody	Negative	Alpha												
CPT: 86022 Patient \$ 160.00				<p><i>Change in Specimen Requirements</i></p> <p>2mL (0.2) Serum Clot Activator SST (Red/Gray or Tiger Top) Refrigerated or Frozen</p> <p>Separate Serum from cells ASAP. Keep Serum refrigerated; however if specimen is not going to be tested within 48 hours then freeze the specimen.</p> <p><i>EDTA Plasma is No Longer Acceptable Specimen Type due to change from reagent vendor.</i></p>																	
Modify Test Effective : 6/5/2017																					
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type												
1501625	D-dimer Quantitative	DDIMR QUAN	Detail	1501625		Ddimr Quan	D-dimer Quantitative	<0.50 ug/mL FEU	Numeric XXX.XX												
CPT: 85379 Patient \$103.70				<p><i>Change in Reference Range and Interpretive Data</i></p> <p>2.7 mL Whole Blood Sodium Citrate 3.2% (Blue Top) Room Temperature or 1.5 mL Double Spun Plasma Sterile, Capped Plastic Tube Frozen</p> <p>Please indicate anticoagulant therapy.</p> <p>Each 2.7mL Sodium Citrate 3.2% (Blue Top) tube must be filled to the proper level, no hemolysis.</p> <p><i>Improperly filled tubes can give erroneous results. Whole blood must be transported to lab immediately.</i></p> <p><i>If testing cannot be started within 4 hours of collection the specimen must be double spun then 1.5 mL plasma aliquot from each tube into individual plastic aliquot tubes and freeze.</i></p> <p>Do not pool aliquots together!</p>																	
				<p><i>Cut off for exclusion of pulmonary embolism (PE) and deep venous thrombosis (DVT) is <0.50 ug/mL FEU</i></p> <p><i>When the quantified D-Dimer value is <0.50 ug/mL FEU and patients have a low or moderate pre-test probability of having venous thromboembolism (PE or DVT), the negative predictive value of this test is 100% for DVT and is 99.7% for PE.</i></p> <p><i>Other tests that can be used in conjunction with the qualitative D-Dimer to diagnose DVT are: compression ultrasound, venography, and MRI. Antithrombotic therapy and prevention of thrombosis, American College of Chest Physician. Evidence-based Clinical Practice Guidelines, 9th edition. 2012</i></p>																	

INACTIVATED TESTS - Please update your EMR catalog with those appropriate to your practice					
INACTIVE Effective : 6/5/2017 - Please update your EMR catalog with those appropriate to your practice					
Test Code	Test Name	Mnemonic	Replacement Information		
5508725	Gelatin Bovine Allergen	GEL BOVINE	Mnemonic	Test Code	Test Name
			GELATIN IC	5613375	Gelatin Allergen (Bovine)
See Build information in the NewTest Section of this Live Wire					
INACTIVE Effective : 6/5/2017 - Please update your EMR catalog with those appropriate to your practice					
Test Code	Test Name	Mnemonic	Replacement Information		
6904891	Borrelia hermsii Antibody Panel, IFA	Borrelia herm Ab	Mnemonic	Test Code	Test Name
			n/a	n/a	n/a
The performing laboratory has discontinued this test and does not offer a replacement assay at this time.					
INACTIVE Effective : 6/30/17 - Please update your EMR catalog with those appropriate to your practice					
Test Code	Test Name	Mnemonic	Replacement Information		
5960200	Trichomonas vaginalis DNA Urine/Swab (Aptima Probe)	TRICH APTI	Mnemonic	Test Code	Test Name
		Urine Specimen	Trich Urine	6910123	Trichomonas Urine
		Swab Specimen	Trich Swab	6910121	Trichomonas Swab
		ThinPrep Specimen	Trich PAP	6910119	Trichomonas PAP
For Replacement Test Build Information See the 05/22/2017 LiveWire.					