

NEW TESTS - Please update your EMR catalog with those appropriate to your practice

The following is a list of new IgE Allergy tests to replace the older Inactivated versions. The new IgG Allergy tests will only report out the numeric count in kU/L
 The collection information, price and CPT codes for each allergen will be the same for each Allergen.

CPT: 86003 Specimen: Separate serum from cells ASAP or within 2 hours of collection.
 Price \$30.05 Transfer 0.25mL(0.1mL) serum per allergen ordered to a Standard Transport Tube.
 Transport Temperature: Refrigerated.
 Unacceptable Conditions: Hemolyzed, icteric, or lipemic specimens.
 Stability After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year
 Set up Sun-Sat Reports 1-2 Days from set up.

New Test Effective : 08/21/2017

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
5613359	Acacia Tree IgE	Acacia Tree	Detail	6015-2	5613359	Acacia Tree IgE	Acacia Tree IgE	< 0.10 kU/L	Numeric XXXX.XX
5613301	Alder Tree IgE	Alder Tree	Detail	15284-3	5613301	Alder Tree IgE	Alder Tree IgE	< 0.10 kU/L	Numeric XXXX.XX
5613303	Alfalfa IgE	Alfalfa	Detail	7065-6	5613303	Alfalfa IgE	Alfalfa IgE	< 0.10 kU/L	Numeric XXXX.XX
5613531	Alpha-Gal (galactose-alpha-1,3-galactose) IgE	Alpha-Gal	Detail	73837-7	5613531	Alpha-Gal IgE	Alpha-Gal IgE	< 0.10 kU/L	Numeric XXXX.XX
5613305	Anchovy IgE	Anchovy	Detail	7079-7	5613305	Anchovy IgE	Anchovy IgE	< 0.10 kU/L	Numeric XXXX.XX
5613307	Asparagus IgE	Asparagus	Detail	7099-5	5613307	Asparagus IgE	Asparagus IgE	< 0.10 kU/L	Numeric XXXX.XX
5613309	Avocado IgE	Avocado	Detail	6033-5	5613309	Avocado IgE	Avocado IgE	< 0.10 kU/L	Numeric XXXX.XX
5613311	Basil IgE	Basil	Detail	7118-3	5613311	Basil IgE	Basil IgE	< 0.10 kU/L	Numeric XXXX.XX
5613313	Bass Black IgE	Black Bass	Detail	7120-9	5613313	Black Bass IgE	Black Bass IgE	< 0.10 kU/L	Numeric XXXX.XX
5613315	Bay Leaf IgE	Bay Leaf	Detail	7125-8	5613315	Bay Leaf IgE	Bay Leaf IgE	< 0.10 kU/L	Numeric XXXX.XX
5613317	Beet Root IgE	Beet Root	Detail	19728-5	5613317	Beet Root IgE	Beet Root IgE	< 0.10 kU/L	Numeric XXXX.XX
5613319	Black Bean IgE	Black Bean	Detail	21096-3	5613319	Black Bean IgE	Black Bean IgE	< 0.10 kU/L	Numeric XXXX.XX
5613479	Black Pepper IgE	Black Pepper	Detail	6045-9	5613479	Black Pepper IgE	Black Pepper IgE	< 0.10 kU/L	Numeric XXXX.XX
5613321	Blackberry IgE	Blackberry	Detail	7151-4	5613321	Blackberry IgE	Blackberry IgE	< 0.10 kU/L	Numeric XXXX.XX
5613323	Blue Mussel IgE	Blue Mussel	Detail	6048-3	5613323	Blue Mussel IgE	Blue Mussel IgE	< 0.10 kU/L	Numeric XXXX.XX
5613325	Blueberry IgE	Blueberry	Detail	7152-2	5613325	Blueberry IgE	Blueberry IgE	< 0.10 kU/L	Numeric XXXX.XX
5613327	Broccoli IgE	Broccoli	Detail	6052-5	5613327	Broccoli IgE	Broccoli IgE	< 0.10 kU/L	Numeric XXXX.XX
5613329	Buckwheat IgE	Buckwheat	Detail	6054-1	5613329	Buckwheat IgE	Buckwheat IgE	< 0.10 kU/L	Numeric XXXX.XX
5613495	Canola(Rapeseed) Oil IgE	Canola Oil	Detail	11195-5	5613495	Canola Oil IgE	Canola Oil IgE	< 0.10 kU/L	Numeric XXXX.XX
5613333	Caraway Seed IgE	Caraway Seed	Detail	7178-7	5613333	Caraway Seed IgE	Caraway Seed IgE	< 0.10 kU/L	Numeric XXXX.XX
5613335	Catfish IgE	Catfish	Detail	7187-8	5613335	Catfish IgE	Catfish IgE	< 0.10 kU/L	Numeric XXXX.XX
5613337	Cauliflower IgE	Cauliflower	Detail	6835-3	5613337	Cauliflower IgE	Cauliflower IgE	< 0.10 kU/L	Numeric XXXX.XX
5613499	Cedar/Red Tree IgE	Red Cedar Tree	Detail	7195-1	5613499	Red Cedar Tr IgE	Red Cedar Tr IgE	< 0.10 kU/L	Numeric XXXX.XX
5613339	Cherry IgE	Cherry	Detail	6719-9	5613339	Cherry IgE	Cherry IgE	< 0.10 kU/L	Numeric XXXX.XX
5613341	Chestnut IgE	Chestnut	Detail	6259-6	5613341	Chestnut IgE	Chestnut IgE	< 0.10 kU/L	Numeric XXXX.XX
5613343	Chickpea (Garbanzo) IgE	Chickpea	Detail	7573-9	5613343	Chickpea IgE	Chickpea IgE	< 0.10 kU/L	Numeric XXXX.XX
5613481	Chili Pepper (C. frutescens) IgE	Chili Pepper	Detail	7591-1	5613481	Chili Pepper IgE	Chili Pepper IgE	< 0.10 kU/L	Numeric XXXX.XX
5613347	Cinnamon IgE	Cinnamon	Detail	6837-9	5613347	Cinnamon IgE	Cinnamon IgE	< 0.10 kU/L	Numeric XXXX.XX
5613349	Cockroach, American IgE	Cockroach, American	Detail	30170-5	5613349	Cockroach Amer IgE	Cockroach Amer IgE	< 0.10 kU/L	Numeric XXXX.XX
5613533	Coffee IgE	Coffee	Detail	6083-0	5613533	Coffee IgE	Coffee IgE	< 0.10 kU/L	Numeric XXXX.XX
5613465	Common Millet IgE	Millet, Common	Detail	6175-4	5613465	Millet IgE	Millet IgE	< 0.10 kU/L	Numeric XXXX.XX
5613345	Coriander/Cilantro IgE	Cilantro/Coriander	Detail	16650-4	5613345	Cilantro/Coriander IgE	Cilantro/Coriander IgE	< 0.10 kU/L	Numeric XXXX.XX
5613351	Cotton Seed IgE	Cotton Seed	Detail	6088-9	5613351	Cotton Seed IgE	Cotton Seed IgE	< 0.10 kU/L	Numeric XXXX.XX
5613353	Cucumber IgE	Cucumber	Detail	6724-9	5613353	Cucumber IgE	Cucumber IgE	< 0.10 kU/L	Numeric XXXX.XX
5613425	Gelatin Porcine IgE	Gelatin Porcine	Detail	16103-4	5613425	Gelatin Porcine IgE	Gelatin Porcine IgE	< 0.10 kU/L	Numeric XXXX.XX
5613427	Ginger IgE	Ginger	Detail	7335-3	5613427	Ginger IgE	Ginger IgE	< 0.10 kU/L	Numeric XXXX.XX
5613429	Goat Milk IgE	Goat Milk	Detail	7340-3	5613429	Goat Milk IgE	Goat Milk IgE	< 0.10 kU/L	Numeric XXXX.XX
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type

5613431	Goldenrod IgE	Goldenrod	Detail	6128-3	5613431	Goldenrod IgE	Goldenrod IgE	< 0.10 kU/L	Numeric XXXX.XX
5613433	Grapefruit IgE	Grapefruit	Detail	6131-7	5613433	Grapefruit IgE	Grapefruit IgE	< 0.10 kU/L	Numeric XXXX.XX
5613435	Green Bean (String) IgE	Green Bean	Detail	6831-2	5613435	Green Bean IgE	Green Bean IgE	< 0.10 kU/L	Numeric XXXX.XX
5613437	Guinea Pig Epithelium IgE	Guinea Pig Epi	Detail	6134-1	5613437	Guinea Pig Epi IgE	Guinea Pig Epi IgE	< 0.10 kU/L	Numeric XXXX.XX
5613423	Gulf Flounder IgE	Flounder, Gulf	Detail	7316-3	5613423	Flounder Gulf IgE	Flounder Gulf IgE	< 0.10 kU/L	Numeric XXXX.XX
5613511	Hackberry IgE	Hackberry	Detail	7383-3	5613511	Hackberry IgE	Hackberry IgE	< 0.10 kU/L	Numeric XXXX.XX
5613439	Haddock IgE	Haddock	Detail	7384-1	5613439	Haddock IgE	Haddock IgE	< 0.10 kU/L	Numeric XXXX.XX
5613441	Halibut IgE	Halibut	Detail	6842-9	5613441	Halibut IgE	Halibut IgE	< 0.10 kU/L	Numeric XXXX.XX
5613443	Hops IgE	Hops	Detail	6141-6	5613443	Hops IgE	Hops IgE	< 0.10 kU/L	Numeric XXXX.XX
5613355	House Dust Greer IgE	House Dust Greer	Detail	9828-5	5613355	HD Greer IgE	HD Greer IgE	< 0.10 kU/L	Numeric XXXX.XX
5613445	Kidney Bean IgE	Kidney Bean	Detail	7129-0	5613445	Kidney Bean IgE	Kidney Bean IgE	< 0.10 kU/L	Numeric XXXX.XX
5613447	Kiwi IgE	Kiwi	Detail	6154-9	5613447	Kiwi IgE	Kiwi IgE	< 0.10 kU/L	Numeric XXXX.XX
5613449	Lemon IgE	Lemon	Detail	6159-8	5613449	Lemon IgE	Lemon IgE	< 0.10 kU/L	Numeric XXXX.XX
5613451	Lentil IgE	Lentil	Detail	7458-3	5613451	Lentil IgE	Lentil IgE	< 0.10 kU/L	Numeric XXXX.XX
5613453	Lime IgE	Lime	Detail	7464-1	5613453	Lime IgE	Lime IgE	< 0.10 kU/L	Numeric XXXX.XX
5613457	Macadamia Nut IgE	Macadamia Nut	Detail	11183-1	5613457	Macadam Nut IgE	Macadamia Nut IgE	< 0.10 kU/L	Numeric XXXX.XX
5613459	Mackerel IgE	Mackerel	Detail	7474-0	5613459	Mackerel IgE	Mackerel IgE	< 0.10 kU/L	Numeric XXXX.XX
5613461	Malt IgE	Malt	Detail	6166-3	5613461	Malt IgE	Malt IgE	< 0.10 kU/L	Numeric XXXX.XX
5613463	Mango IgE	Mango	Detail	6167-1	5613463	Mango IgE	Mango IgE	< 0.10 kU/L	Numeric XXXX.XX
5613467	Mouse Urine IgE	Mouse Urine	Detail	6181-2	5613467	Mouse Ur IgE	Mouse Ur IgE	< 0.10 kU/L	Numeric XXXX.XX
5613469	Mushroom IgE	Mushroom	Detail	6848-6	5613469	Mushroom IgE	Mushroom IgE	< 0.10 kU/L	Numeric XXXX.XX
5613471	Mustard IgE	Mustard	Detail	6185-3	5613471	Mustard IgE	Mustard IgE	< 0.10 kU/L	Numeric XXXX.XX
5613487	Navy Bean IgE	Navy Bean	Detail	6279-4	5613487	Navy Bean IgE	Navy Bean IgE	< 0.10 kU/L	Numeric XXXX.XX
5613473	Nutmeg IgE	Nutmeg	Detail	6188-7	5613473	Nutmeg IgE	Nutmeg IgE	< 0.10 kU/L	Numeric XXXX.XX
5613475	Orchard Grass/Cocksfoot IgE	Orchard Grass	Detail	6195-2	5613475	Orchard Grass IgE	Orchard Grass IgE	< 0.10 kU/L	Numeric XXXX.XX
5613477	Oregano IgE	Oregano	Detail	7552-3	5613477	Oregano IgE	Oregano IgE	< 0.10 kU/L	Numeric XXXX.XX
5613535	Pepper C. annuum IgE	Pepper C Annum	Detail	7586-1	5613535	Pepper C annum IgE	Pepper C annum IgE	< 0.10 kU/L	Numeric XXXX.XX
5613483	Perch IgE	Perch	Detail	15930-1	5613483	Perch IgE	Perch IgE	< 0.10 kU/L	Numeric XXXX.XX
5613485	Pine (Pinon) Nut IgE	Pine Nut	Detail	10954-6	5613485	Pine Nut IgE	Pine Nut IgE	< 0.10 kU/L	Numeric XXXX.XX
5613489	Plum IgE	Plum	Detail	6853-6	5613489	Plum IgE	Plum IgE	< 0.10 kU/L	Numeric XXXX.XX
5613493	Rabbit Epithelium IgE	Rabbit Epithelium	Detail	6223-2	5613493	Rabbit Epi IgE	Rabbit Epi IgE	< 0.10 kU/L	Numeric XXXX.XX
5613497	Raspberry IgE	Raspberry	Detail	7656-2	5613497	Raspberry IgE	Raspberry IgE	< 0.10 kU/L	Numeric XXXX.XX
5613357	Red Dye (Carmine) IgE	Carmine	Detail	17511-7	5613357	Carmine IgE	Carmine IgE	< 0.10 kU/L	Numeric XXXX.XX
5613501	Sardine IgE	Sardine	Detail	11197-1	5613501	Sardine IgE	Sardine IgE	< 0.10 kU/L	Numeric XXXX.XX
5613503	Spinach IgE	Spinach	Detail	6856-9	5613503	Spinach IgE	Spinach IgE	< 0.10 kU/L	Numeric XXXX.XX
5613507	Squid IgE	Squid	Detail	14038-4	5613507	Squid IgE	Squid IgE	< 0.10 kU/L	Numeric XXXX.XX
5613331	Sugar Cane IgE	Sugar Cane	Detail	7723-0	5613331	Sugar Cane IgE	Sugar Cane IgE	< 0.10 kU/L	Numeric XXXX.XX
5613505	Summer Squash IgE	Squash, Summer	Detail	7714-9	5613505	Squash IgE	Squash IgE	< 0.10 kU/L	Numeric XXXX.XX
5613529	Sweet Gum Tree IgE	Sweet Gum Tree	Detail	6260-4	5613529	Sweet Gum Tr IgE	Sweet Gum Tr IgE	< 0.10 kU/L	Numeric XXXX.XX
5613509	Tea IgE	Tea	Detail	6264-6	5613509	Tea IgE	Tea IgE	< 0.10 kU/L	Numeric XXXX.XX
5613513	Trout IgE	Trout	Detail	6268-7	5613513	Trout IgE	Trout IgE	< 0.10 kU/L	Numeric XXXX.XX
5613515	Vanilla IgE	Vanilla	Detail	7756-0	5613515	Vanilla IgE	Vanilla IgE	< 0.10 kU/L	Numeric XXXX.XX
5613517	Whey IgE	Whey	Detail	7774-3	5613517	Whey IgE	Whey IgE	< 0.10 kU/L	Numeric XXXX.XX
5613519	White Pine Tree IgE	White Pine Tree	Detail	6282-8	5613519	White Pine Tr IgE	White Pine Tr IgE	< 0.10 kU/L	Numeric XXXX.XX
5613521	Whitefish IgE	Whitefish	Detail	7775-0	5613521	Whitefish IgE	Whitefish IgE	< 0.10 kU/L	Numeric XXXX.XX
5613523	Willow Tree IgE	Willow Tree	Detail	6285-1	5613523	Willow Tr IgE	Willow Tr IgE	< 0.10 kU/L	Numeric XXXX.XX
5613525	Yellow Dock IgE	Yellow Dock	Detail	6738-9	5613525	Yellow Dock IgE	Yellow Dock IgE	< 0.10 kU/L	Numeric XXXX.XX
5613527	Zucchini IgE	Zucchini	Detail	7718-0	5613527	Zucchini IgE	Zucchini IgE	< 0.10 kU/L	Numeric XXXX.XX

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
5613551	Allergen Interp, Immucap Score IgE	Allergen Score IgE	Detail	33536-4	5613551	Allergen Score IgE	Allergen Score IgE	By Report	Alpha

This test will be added on to each encounter Score. to report the interpretive data to help you Crossreference the numerical value of kU/L to the corresponding Allergy Class

Allergen, Interp, Immucap Score IgE:

Less than 0.10 kU/L	Class 0	No significant level detected
0.10-0.34 kU/L	Class 0/1	Clinical relevance undetermined
0.35-0.70 kU/L	Class 1	Low
0.71-3.50 kU/L	Class 2	Moderate
3.51-17.50 kU/L	Class 3	High
17.51-50.00 kU/L	Class 4	Very High
50.01-100.00 kU/L	Class 5	Very High

MODIFICATION TESTS - Please update your EMR catalog with those appropriate to your practice
Modify Test Effective : 08/21/2017
- Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	Careset Component Information						
2929745	Weed Allergy Panel (Custom)	AO WEED PN	Care Set	Order Code	Result Code	Mnemonic	Test Name	Reference Range UOM	Result Type	
CPT: 86003x8 Price: \$240.40 2.0mL (0.8mL) Serum Clot Activator SST (Red/Gray or Tiger Top) Specimen: Separate serum from cells ASAP or within 2 hours of collection. Transport/Storage: Refrigerated Specimen Stability: refrigerated 4 weeks, frozen 4 weeks. Unacceptable: Room Temperature 5606925 Allergen Interp * Only one of these Interpretations will be reported per panel				5644025		PIGWEEED	Common Pigweed Allergen			
					5644075	Pigweed Cl	Pigweed Class	Class 0	Alpha	
					5644050	Pigweed Ct	Pigweed Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX	
					5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha	
					2930700	YELLOWDOG	Yellowdock Allergen			
					2930650	Yellowdock Cl	Yellowdock Class			
					5533600	Yellowdock Ct	Yellowdock Count kU/L			
					5613525	5613525	Yellow Dock	Yellow Dock IgE	<0.10 kU/L	Numeric XXXX.XX
					2930500	RAGWEED W	Western Ragweed Allergen			
					3630450	West Ragweed Cl	West Ragweed Class	Class 0	Alpha	
					5527900	West Ragweed Ct	West Ragweed Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX	
					5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha	
					5609225	COCKLEBUR	Cocklebur Allergen			
					5609200	Cocklebur Cl	Cocklebur Class	Class 0	Alpha	
					5609175	Cocklebur Ct	Cocklebur Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX	
					5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha	
					5609300	LAMBSQUART	Goosefoot, Lamb's Quarters Allergen			
					5609275	Goosefoot Cl	Goosefoot Class	Class 0	Alpha	
					5609250	Goosefoot Ct	Goosefoot Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX	
					5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha	
					5650125	NETTLE	Nettle Allergen			
					5650075	Nettle Cl	Nettle Class	Class 0	Alpha	
					5650025	Nettle Ct	Nettle Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX	
					5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha	
					5607625	ENGL PLANT	Plantain (English), Ribwort Allergen			
	5607600	Plantain Rib Cl	Plantain Rib Class	Class 0	Alpha					
	5607575	Plantain Rib Ct	Plantain Rib Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX					
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha					
	5612025	MUGWORT IC	Mugwort Allergen							
	5612000	Mugwort Cl	Mugwort Class	Class 0	Alpha					
	5611975	Mugwort Ct	Mugwort Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX					
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha					

Modify Test Effective : 08/21/2017
- Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	Careset Component Information							
3492170	Allergen Profile Zone 7	A ZONE 7	Care Set	Component	Result Code	Mnemonic	Test Name	Reference Range	UOM	Result Type	
CPT: 86003x31 Price \$931.55 10mL (3mL) Serum Clot Activator SST (Red/Gray or Tiger Top) Specimen: Separate serum from cells ASAP or within 2 hours of collection. Transport/Storage: Refrigerated Specimen Stability: refrigerated 4 weeks, frozen 4 weeks. Unacceptable: Room Temperature 5606925 Allergen Interp * Only one of these Interpretations will be reported per panel				5606550		ALTER TEN	Alternaria alternata Allergen				
					5606525	Alt Ten Cl	Alt Ten Class	Class 0		Alpha	
					5606500	Alt Ten Ct	Alt Ten Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX	
						5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha
					5606400		ASPER MLD	Aspergillus fumigatus Allergen			
						5606375	Asp Fum Cl	Asp Fum Class	Class 0		Alpha
						5606350	Asp Fum Ct	Asp Fum Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX
						5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha
					5617075		BAHIA GRS	Bahia Grass Allergen			
						5617050	Bahia Grass Cl	Bahia Grass Class	Class 0		Alpha
						5617025	Bahia Grass Ct	Bahia Grass Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX
						5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha
					5606325		BERMUDA	Bermuda Grass Allergen			
						5606300	Bermuda Cl	Bermuda Class	Class 0		Alpha
						5606275	Bermuda Ct	Bermuda Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX
						5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha
					5611750		BOX ELDER	Box Elder (Maple) Tree Allergen			
						5611725	Box-elder Cl	Box-elder Class	Class 0		Alpha
						5611700	Box-elder Ct	Box-elder Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX
						5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha
					5606250		CAT DANDER	Cat Dander Allergen			
						5606225	Cat Dander Cl	Cat Dander Class	Class 0		Alpha
						5606200	Cat Dander Ct	Cat Dander Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX
						5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha
					5606175		CLADOSPOR	Classadosporium herbarum Allergen			
						5606150	Cladosporium Cl	Classadosporium Class	Class 0		Alpha
						5606125	Cladosporium Ct	Classadosporium Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX
						5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha
					5605800		DOG DAND	Dog Dander Allergen			
						5605750	Dog Dander Cl	Dog Dander Class	Class 0		Alpha
						5605775	Dog Dander Ct	Dog Dander Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX
						5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha
					5604985		MITE DUST	Dust Mite (D. farinae) Allergen			
					5604975	Dust Mite Cl	Dust Mite Class	Class 0		Alpha	
					5604950	Dust Mite Ct	Dust Mite Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX	
					5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha	
				5606795		ELM TREE	Elm Tree Allergen				
					5506785	Elm Tree Cl	Elm Tree Class	Class 0		Alpha	
					5606775	Elm Tree Ct	Elm Tree Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX	
					5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha	
				5606025		COCKROACH	Cockroach, German Allergen				
					5606000	Cockroach Cl	Cockroach Class	Class 0		Alpha	
					5605975	Cockroach Ct	Cockroach Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX	
					5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha	

Allergen Profile Zone 7 Continued Next Page

Allergen Profile Zone 7 Continued

5609675		HOUSE DUST	House Dust Allergen		
	5609650	House Dust Cl	House Dust Class	Class 0	Alpha
	5609625	House Dust Ct	House Dust Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5607550		DERM PTE	House Dust Mite Allergen		
	5607525	Derm pterony Cl	Derm pterony Class	Class 0	Alpha
	5607500	Derm pterony Ct	Derm pterony Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5527975		JOHNS SMUT	Johnson Smut Grass Allergen		
	5527775	John Smut Gr Cl	John Smut Gr Class	Class 0	Alpha
	5527925	John Smut Gr Ct	John Smut Gr Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
5605885		KENT BLUE	Meadow Grass, Kentucky Blue		
	5605875	Kent Blue Cl	Kent Blue Class	Class 0	Alpha
	5605850	Kent Blue Ct	Kent Blue Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5613225		MESQUITE	Mesquite Tree Allergen		
	5613200	Mesquite Tr Cl	Mesquite Tr Class	Class 0	Alpha
	5613175	Mesquite Tr Ct	Mesquite Tr Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5611850		WE JUNIPR	Mountain Juniper Tree Allergen		
	5611825	Mtn Juniper Cl	Mtn Juniper Class	Class 0	Alpha
	5611800	Mtn Juniper Ct	Mtn Juniper Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5621325		MUCOR RACE	Mucor racemosus Allergen		
	5621300	Mucor Race Cl	Mucor Race Class	Class 0	Alpha
	5621275	Mucor Race Ct	Mucor Race Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5518175		MULBERRY RED	Mullberry Tree Red Allergen		
	5518150	Mulberry Red Cl	Mulberry Red Class	Class 0	Alpha
	5528650	Mulberry Red Ct	Mulberry Red Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5650125		NETTLE	Nettle Allergen		
	5650075	Nettle Cl	Nettle Class	Class 0	Alpha
	5650025	Nettle Ct	Nettle Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5604925		OAK TREE	Oak Tree Allergen		
	5590800	Oak Tree Cl	Oak Tree Class	Class 0	Alpha
	5528800	Oak Tree Ct	Oak Tree Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5611150		PECAN MEAT	Pecan Allergen		
	5611125	Pecan Meat Cl	Pecan Meat Class	Class 0	Alpha
	5611100	Pecan Meat Ct	Pecan Meat Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5622400		PEN CHRYSO	Penicillin Chrysogonium Allergen		
	5622375	Penicillium Cl	Penicillium Class	Class 0	Alpha
	5622350	Penicillium Ct	Penicillium Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5644025		PIGWEED	Common Pigweed Allergen		
	5644075	Pigweed Cl	Pigweed Class	Class 0	Alpha
	5644050	Pigweed Ct	Pigweed Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha

Allergen Profile Zone 7 Continued Next Page

Allergen Profile Zone 7 Continued

5607625		ENGL PLANT	Plantain (English), Ribwort Allergen		
	5607600	Plantain Rib Cl	Plantain Rib Class	Class 0	Alpha
	5607575	Plantain Rib Ct	Plantain Rib Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5609075		SHEEP SORL	Sheep Sorrel Allergen		
	5609050	Sheep Sorrel Cl	Sheep Sorrel Class	Class 0	Alpha
	5609025	Sheep Sorrel Ct	Sheep Sorrel Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5621625		STEMPHYLLI	Stemphylium botryosum Allergen		
	5621600	Stemphylium Cl	Stemphylium Class	Class 0	Alpha
	5621575	Stemphylium Ct	Stemphylium Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
	5518840	SWEETGUM	Sweet Gum Tree Allergen		
	5518830	Swt Gum Tr Cl	Swt Gum Tr Class		
	5518835	Swt Gum Tr Ct	Swt Gum Tr Count kU/L		
	5613529	Sweet Gum Tree	Sweet Gum Tree IgE	<0.10 kU/L	Numeric XXXX.XX
2930500		RAGWEED W	Western Ragweed Allergen		
	3630450	West Ragweed Cl	West Ragweed Class	Class 0	Alpha
	5527900	West Ragweed Ct	West Ragweed Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5606475		ASH WHITE	White Ash Tree Allergen		
	5606450	Wht Ash Tr Cl	Wht Ash Tr Class	Class 0	Alpha
	5606425	Wht Ash Tr Ct	Wht Ash Tr Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5518075		WINGSCALE	Wingscale, Lenscale Allergen		
	5518050	Wingscale Cl	Wingscale Class	Class 0	Alpha
	5528000	Wingscale Ct	Wingscale Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha

Remove This Code

Add This Replacement Code

Modify Test Effective : 08/21/2017

- Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	Careset Component Information						
2929945	Insect Allergy Panel (Custom)	AO INSECT	Care Set	Component	Result Code	Mnemonic	Test Name	Reference Range UOM	Result Type	
CPT: 86003x4 Price \$120.20 1.0mL (0.5mL) Serum Clot Activator SST (Red/Gray or Tiger Top) Specimen: Separate serum from cells ASAP or within 2 hours of collection. Transport/Storage: Refrigerated Specimen Stability: refrigerated 4 weeks, frozen 4 weeks. Unacceptable: Room Temperature 5606925 Allergen Interp * Only one of these Interpretations will be reported per panel				5607550		DERM PTE	House Dust Mite Allergen			
					5607525	Derm pterony Cl	Derm pterony Class	Class 0	Alpha	
					5607500	Derm pterony Ct	Derm pterony Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX	
					5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha	
					5606025		COCKROACH	Cockroach, German Allergen		
					5606000	Gr Roach CL	Gr Roach CL	Class 0	Alpha	
					5605975	Gr Roach CT	Gr Roach CT	<0.3.5 kU/L	Numeric XXXX.XX	
					5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha	
					5616885		COCKRO-AM	Cockroach, American Allergen		
					5618645		Amer-Roach Ct	Amer-Roach Count kU/L		
					5616945		Amer-Roach Cl	Amer-Roach Class		
					5613349	5613349	Cockroach, American	Cockroach, American IgE	<0.10 kU/L	Numeric XXXX.XX
					5604985		MITE DUST	Dust Mite (D. farinae) Allergen		
					5604975	Dust Mite Cl	Dust Mite Class	Class 0	Alpha	
					5604950	Dust Mite Ct	Dust Mite Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX	
				5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha		

Remove This Code

Add This Replacement Code

Modify Test Effective : 08/21/2017

- Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	Careset Component Information							
2929545	Tree Allergy Panel (Custom)	AO TREE PN	Care Set	Component	Result Code	Mnemonic	Test Name	Reference Range	UOM	Result Type	
CPT: 86003x14 Price \$420.70 2.0mL (1.5mL) Serum Clot Activator SST (Red/Gray or Tiger Top) Specimen: Separate serum from cells ASAP or within 2 hours of collection. Transport/Storage: Refrigerated Specimen Stability: refrigerated 4 weeks, frozen 4 weeks. Unacceptable: Room Temperature 5606925 Allergen Interp * Only one of these Interpretations will be reported per panel				5606475		ASH WHITE	White Ash Tree Allergen				
					5606450	Wht Ash Tr Cl	Wht Ash Tr Class	Class 0		Alpha	
					5606425	Wht Ash Tr Ct	Wht Ash Tr Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX	
					5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha	
					5611750		BOX ELDER	Box Elder (Maple) Tree Allergen			
						5611725	Box-elder Cl	Box-elder Class	Class 0		Alpha
						5611700	Box-elder Ct	Box-elder Count kU/L	<0.3.5 kU/L		Numeric XXXX.XX
						5606925	Allergen Interp *	Allergen Interpation *	By Report		Alpha
					2954175		RED CEDAR	Cedar Tree Allergen			
						2954150	Cedar Tree Cl	Cedar Tree Class			
						2954125	Cedar Tree Ct	Cedar Tree Count kU/L			
					Add This Replacement Code	5613499	5613499	Red Cedar Tree	Cedar/Red Tree IgE	<0.10 kU/L	Numeric XXXX.XX
						5607475		COTONWOOD	Cottonwood Tree Allergen		
							5607450	Cottonwood Cl	Cottonwood Class	Class 0	Alpha
							5607425	Cottonwood Ct	Cottonwood Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
							5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
					Remove This Code	5518525		WILLOW TR	Willow Pine Tree Allergen		
							5518500	Willow Pine Cl	Willow Pine Class		
							5528350	Willow Pine Ct	Willow Pine Count kU/L		
					Add This Replacement Code	5613523	5613523	Willow Tree	Willow Tree IgE	<0.10 kU/L	Numeric XXXX.XX
					Remove This Code	5520850		TREE HACK	Hackberry Tree Allergen		
							5518100	Hackberry Cl	Hackberry Tree Allergen		
							5528600	Hackberry Ct	Hackberry Tree Allergen		
					Add This Replacement Code	5613511	5613511	Hackberry	Hackberry IgE	<0.10 kU/L	Numeric XXXX.XX
						5606795		ELM TREE	Elm Tree Allergen		
							5506785	Elm Tree Cl	Elm Tree Class	Class 0	Alpha
							5606775	Elm Tree Ct	Elm Tree Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
							5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
						5606750		PECAN TREE	Pecan, Hickory Tree Allergen		
							5606725	Pecan Tr Cl	Pecan Tr Class	Class 0	Alpha
						5606700	Pecan Tr Ct	Pecan Tr Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX	
						5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha	
					5604475		WALNUT TRE	Walnut Tree Allergen			
						5604450	Walnut Tr Cl	Walnut Tr Class	Class 0	Alpha	
						5604425	Walnut Tr Ct	Walnut Tr Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX	
						5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha	
					5518175		MULBERRY RED	Mullberry Tree Red Allergen			
						5518150	Mulberry Red Cl	Mulberry Red Class	Class 0	Alpha	
						5528650	Mulberry Red Ct	Mulberry Red Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX	
						5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha	
					5604925		OAK TREE	Oak Tree Allergen			
						5590800	Oak Tree Cl	Oak Tree Class	Class 0	Alpha	
						5528800	Oak Tree Ct	Oak Tree Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX	
						5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha	

Tree Allergy Panel (Custom) Continued Next Page

Tree Allergy Panel (Custom) Continued

5613300		SYCAMRE IC	Sycamore, London Plane Tree Allergen		
	5613275	Sycamore CI	Sycamore Class	Class 0	Alpha
	5613250	Sycamore Ct	Sycamore Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5613000		SILV BIRCH	Common Silver Birch Tree Allergen		
	5612075	Silver Birch CI	Silver Birch Class	Class 0	Alpha
	5612050	Silver Birch Ct	Silver Birch Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha
5613225		MESQUITE	Mesquite Tree Allergen		
	5613200	Mesquite Tr CI	Mesquite Tr Class	Class 0	Alpha
	5613175	Mesquite Tr Ct	Mesquite Tr Count kU/L	<0.3.5 kU/L	Numeric XXXX.XX
	5606925	Allergen Interp *	Allergen Interpation *	By Report	Alpha

INACTIVATED TESTS - Please update your EMR catalog with those appropriate to your practice

This section is a cross reference list of Allergy tests to be Inactivated associated to the New IgE Allergy tests for replacement
The Replacement test build information is listed at the top Green Section of the LiveWire.

INACTIVE Effective : 08/21/2017		Replacement Information				
Inactive	Inactive	Inactive		Replacement	Replacement	Replacement
Test Code	Test Name	Mnemonic		Test Code	Mnemonic	Test Name
5594260	Acacia Tree Allergen	ACASIA	Replacement	Test Code: 5613359	Acacia Tree	Acacia Tree IgE
5518225	Alder Red (Alnus Rubra) Allergen	ALDR T TR	Replacement	Test Code: 5613301	Alder Tree	Alder Tree IgE
5560825	Alfalfa Allergen	ALFALFA	Replacement	Test Code: 5613303	Alfalfa	Alfalfa IgE
5519675	Alpha 1,3 Galactose IgE	ALPHA GAL	Replacement	Test Code: 5613531	Alpha-Gal	Alpha-Gal (galactose-alpha-1,3-galactose) IgE
5552350	Anchovy Allergen	ANCHOVY	Replacement	Test Code: 5613305	Anchovy	Anchovy IgE
5556025	Asparagus Allergen	ASPARAGUS	Replacement	Test Code: 5613307	Asparagus	Asparagus IgE
5595050	Avocado Allergen	AVOCADO	Replacement	Test Code: 5613309	Avocado	Avocado IgE
5517175	Basil Allergen	BASIL	Replacement	Test Code: 5613311	Basil	Basil IgE
5518125	Bass Allergen	BASS	Replacement	Test Code: 5613313	Black Bass	Bass Black IgE
5516190	Bayleaf allergen	BAYLEAF	Replacement	Test Code: 5613315	Bay Leaf	Bay Leaf IgE
5555875	Beets (Beetroot) Allergen	BEET	Replacement	Test Code: 5613317	Beet Root	Beet Root IgE
5559735	Black Bean Allergen	BLACK BEAN	Replacement	Test Code: 5613319	Black Bean	Black Bean IgE
5532450	Black Pepper Allergen	PEPPER BLK	Replacement	Test Code: 5613479	Black Pepper	Black Pepper IgE
5561375	Blackberry Allergen	BLACKBERRY	Replacement	Test Code: 5613321	Blackberry	Blackberry IgE
5566675	Blue Mussel Allergen	BLUE MUSSL	Replacement	Test Code: 5613323	Blue Mussel	Blue Mussel IgE
5594000	Blueberry Allergen	BLUEBERRY	Replacement	Test Code: 5613325	Blueberry	Blueberry IgE
5597000	Broccoli Allergen	BROCCOLI	Replacement	Test Code: 5613327	Broccoli	Broccoli IgE
5556500	Buckwheat Allergen	BUCKWHEAT	Replacement	Test Code: 5613329	Buckwheat	Buckwheat IgE
5536875	Rape Seed (Canola) Allergen	RAPE SEED	Replacement	Test Code: 5613495	Canola Oil	Canola(Rapeseed) Oil IgE
5500335	Caraway Allergen	CARAWAY SD	Replacement	Test Code: 5613333	Caraway Seed	Caraway Seed IgE
5610475	Catfish Allergen	CATFISH	Replacement	Test Code: 5613335	Catfish	Catfish IgE
5506290	Cauliflower Allergen	CAULIFLOWR	Replacement	Test Code: 5613337	Cauliflower	Cauliflower IgE
2954175	Cedar Tree Allergen	RED CEDAR	Replacement	Test Code: 5613499	Red Cedar Tree	Cedar/Red Tree IgE
5522200	Cherry Allergen	CHERRY	Replacement	Test Code: 5613339	Cherry	Cherry IgE
5500531	Chestnut Food Allergen	CHESTNT FD	Replacement	Test Code: 5613341	Chestnut	Chestnut IgE
5559175	Chickpea Allergen	CHICK PEA	Replacement	Test Code: 5613343	Chickpea	Chickpea (Garbanzo) IgE
5533450	Pepper, Chili Allergen	PEPPR CHIL	Replacement	Test Code: 5613481	Chili Pepper	Chili Pepper (C. frutescens) IgE
5592865	Cinnamon Allergen	CINNAMON	Replacement	Test Code: 5613347	Cinnamon	Cinnamon IgE
~ Continued on Next Page ~						
Inactive	Inactive	Inactive		Replacement	Replacement	Replacement

Test Code	Test Name	Mnemonic		Test Code	Mnemonic	Test Name
5616885	Cockroach, American Allergen	COCKRO AM	Replacement	5613349	Cockroach, American	Cockroach, American IgE
5650135	Coffee Allergen	COFFEE IGE	Replacement	5613533	Coffee	Coffee IgE
5506160	Common Millet Allergen	MILLET	Replacement	5613465	Millet, Common	Common Millet IgE
5516175	Cilantro and Coriander Allergen	CILANTRO	Replacement	5613345	Cilantro/Coriander	Coriander/Cilantro IgE
5556400	Cotton Seed Allergen	COTTNSEED	Replacement	5613351	Cotton Seed	Cotton Seed IgE
5578400	Cucumber Allergen	CUCUMBER	Replacement	5613353	Cucumber	Cucumber IgE
5508575	Gelatin Porcine (Pig) Allergen	GEL PORCIN	Replacement	5613425	Gelatin Porcine	Gelatin Porcine IgE
5558435	Ginger Allergen	GINGER	Replacement	5613427	Ginger	Ginger IgE
6900709	Goat Milk Allergen IgE	Goat Milk IgE	Replacement	5613429	Goat Milk	Goat Milk IgE
5517925	Goldenrod Allergen	GOLDENROD	Replacement	5613431	Goldenrod	Goldenrod IgE
5557825	Grapefruit Allergen	GRAPEFRUIT	Replacement	5613433	Grapefruit	Grapefruit IgE
5597075	Green Bean Allergen	GREEN BEAN	Replacement	5613435	Green Bean	Green Bean (String) IgE
3660250	Guinea Pig Epithelium Allergen	GUINEA PIG	Replacement	5613437	Guinea Pig Epi	Guinea Pig Epithelium IgE
5599075	Flounder Allergen	FLOUNDER	Replacement	5613423	Flounder, Gulf	Gulf Flounder IgE
5520850	Hackberry Tree Allergen	TREE HACK	Replacement	5613511	Hackberry	Hackberry IgE
5599185	Haddock Allergen	HADDOCK	Replacement	5613439	Haddock	Haddock IgE
5552250	Halibut Allergen	HALIBUT	Replacement	5613441	Halibut	Halibut IgE
5519825	Hops Allergen	HOPS	Replacement	5613443	Hops	Hops IgE
5579650	Dust Greer Allergen	DUST GREER	Replacement	5613355	House Dust Greer	House Dust Greer IgE
5515325	Red Kidney Bean Allergen	KIDNY BEAN	Replacement	5613445	Kidney Bean	Kidney Bean IgE
5532075	Kiwi Allergen	KIWI	Replacement	5613447	Kiwi	Kiwi IgE
5594075	Lemon Allergen	LEMON	Replacement	5613449	Lemon	Lemon IgE
5506230	Lentil Allergen	LENTIL	Replacement	5613451	Lentil	Lentil IgE
5518825	Lime Allergen	LIME	Replacement	5613453	Lime	Lime IgE
5500425	Macadamia Nut Allergen	MACADAMIA	Replacement	5613457	Macadamia Nut	Macadamia Nut IgE
6900717	Mackerel Allergen IgE	Mackerel Allergen	Replacement	5613459	Mackerel	Mackerel IgE
5620175	Malt Allergen	MALT	Replacement	5613461	Malt	Malt IgE
5519035	Mango Allergen	MANGO	Replacement	5613463	Mango	Mango IgE
5607090	Mouse Urine Proteins Allergen	MOUSE URIN	Replacement	5613467	Mouse Urine	Mouse Urine IgE
5598090	Mushroom Allergen	MUSHROOM	Replacement	5613469	Mushroom	Mushroom IgE
5562750	Mustard Seed Allergen	MUSTARD	Replacement	5613471	Mustard	Mustard IgE
5582750	White Bean Allergen	WHITE BEAN	Replacement	5613487	Navy Bean	Navy Bean IgE
5598050	Pinto Bean Allergen	PINTO BEAN	Replacement	5613487	Navy Bean	Navy Bean IgE
5515835	Nutmeg Allergen	NUTMEG	Replacement	5613473	Nutmeg	Nutmeg IgE
3633050	Orchard Grass Allergen	ORCHRD GR	Replacement	5613475	Orchard Grass	Orchard Grass/Cocksfoot IgE
5516875	Oregano Allergen	OREGANO	Replacement	5613477	Oregano	Oregano IgE
5532900	Cayenne Pepper Allergen	PEPPER CAY	Replacement	5613535	Pepper C Annum	Pepper C. annum IgE
5536125	Jalapeno Pepper Allergen	PEPPER JAL	Replacement	5613535	Pepper C Annum	Pepper C. annum IgE
5578700	Pepper, Green Allergen	BELL PEPPE	Replacement	5613535	Pepper C Annum	Pepper C. annum IgE
5536275	Perch Allergen	PERCH	Replacement	5613483	Perch	Perch IgE
5552220	Pinion Nut Allergen	PINION NUT	Replacement	5613485	Pine Nut	Pine (Pinon) Nut IgE
5506110	Plum Allergen	PLUM	Replacement	5613489	Plum	Plum IgE
5533850	Rabbit Epithelium Allergen	RABBIT EPI	Replacement	5613493	Rabbit Epithelium	Rabbit Epithelium IgE
5506370	Raspberry Allergen	RASBERRY	Replacement	5613497	Raspberry	Raspberry IgE
5525125	Carmine Allergen	CARMINE	Replacement	5613357	Carmine	Red Dye (Carmine) IgE
5556325	Sardine Allergen	SARDINE	Replacement	5613501	Sardine	Sardine IgE
5506400	Spinach Allergen	SPINACH	Replacement	5613503	Spinach	Spinach IgE

~ Continued on Next Page ~

Inactive	Inactive	Inactive			Replacement	Replacement	Replacement
Test Code	Test Name	Mnemonic			Test Code	Mnemonic	Test Name
5594297	Pacific Squid Allergen	SQUID PACF	Replacement	Test Code:	5613507	Squid	Squid IgE
5594291	Squid Allergen	SQUID IGE	Replacement	Test Code:	5613507	Squid	Squid IgE
5515475	Cane Sugar Allergen	CANE SUGAR	Replacement	Test Code:	5613331	Sugar Cane	Sugar Cane IgE
5555725	Squash Allergen	SQUASH	Replacement	Test Code:	5613505	Squash, Summer	Summer Squash IgE
5518840	Sweet Gum Tree Allergen	SWEETGUM	Replacement	Test Code:	5613529	Sweet Gum Tree	Sweet Gum Tree IgE
6900727	Tea Allergen IgE	Tea IgE	Replacement	Test Code:	5613509	Tea	Tea IgE
5552200	Trout Allergen	TROUT	Replacement	Test Code:	5613513	Trout	Trout IgE
5533625	Vanilla Allergen	VANILLA	Replacement	Test Code:	5613515	Vanilla	Vanilla IgE
5583400	Whey Allergen	WHEY	Replacement	Test Code:	5613517	Whey	Whey IgE
5515450	White Pine Tree Allergen	WHITE PIN	Replacement	Test Code:	5613519	White Pine Tree	White Pine Tree IgE
5536625	Whitefish Allergen	WHITEFISH	Replacement	Test Code:	5613521	Whitefish	Whitefish IgE
5518525	Willow Pine Tree Allergen	WILLOW TR	Replacement	Test Code:	5613523	Willow Tree	Willow Tree IgE
2930700	Yellowdock Allergen	YELLOWDOC	Replacement	Test Code:	5613525	Yellow Dock	Yellow Dock IgE
5556625	Zucchini Allergen	ZUCCHINI	Replacement	Test Code:	5613527	Zucchini	Zucchini IgE