

NEW TESTS - Please update your EMR catalog with those appropriate to your practice

New Test Effective : 11/27/2017					
Test Code	Test Name	Mnemonic	Category/Type	Careset Component Information	
4000334	Pain Panel Comprehensive Urine	PAIN PNLC	Care Set		
CPT: See Components Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			Order Code	Mnemonic	Test Name
			4000107	.SC ADULTR	.Adulterants Screen
			4000100	.Prescriptions	.Prescriptions
			4000335	.PAIN PNLC	.Pain Panel Comprehensive

New Test Effective : 11/27/2017					
Test Code	Test Name	Mnemonic	Category/Type	Careset Component Information	
4000314	Pain Panel General Urine	PAIN PNLG	Care Set		
CPT: See Components Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			Order Code	Mnemonic	Test Name
			4000107	.SC ADULTR	.Adulterants Screen
			4000100	.Prescriptions	.Prescriptions
			4000315	.PAIN PNLG	.Pain Panel General

New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000107	.Adulterants Screen Urine	.SC ADULTR	Group						
CPT: N/A Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months This test should only be ordered as a component of the Pain Panel Comprehensive Urine (4000334) or Pain Panel General Urine (4000314)				2161-8	4000074	Creatinine U	Creatinine Urine	20 - 400 mg/dL	XXXX
				2756-5	4000072	pH U	pH Urine	4.00 - 10.00	XX.XX
				2965-2	4000073	Specific Grav U	Specific Gravity Urine	1.003 - 1.025	X.XXX
				N/A	4000071	Adulterants	Adulterants Screen	N/A	Alpha

New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000106	.Prescriptions	.Prescriptions	Detail	N/A	4000106	Prescriptions	Prescriptions	N/A	Textual
CPT: N/A This is an information test required as part of the other orders on this document, there is no specimen required for this test. This test should only be ordered as a component of the Pain Panel Comprehensive Urine (4000334) or Pain Panel General Urine (4000314)									

New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000335	.Pain Panel Comprehensive Urine	.PAIN PNLC	Group						
CPT: 80307 HCPCS: G0479 Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months A "Detected" result on the following screen components will reflex the indicated confirmation order: Screen (assay): (400070)Acetaminophen (400051)Amphetamines (400052)Barbiturates (400053)Benzodiazepines (400060)Buprenorphine (400081)Cannabinoids (400056)Carisoprodol (400080)Cocaine (400062)Fentanyl (400059)Gabapentin (400064)Methadone (400084)Methylphenidate (400054)Opiates (400079)Phencyclidine (400057)Pregabalin (400066)Propoxyphene (400083)Pseudoephedrine (400068)Tapentadol (400048)Tramadol (400082)Trazodone (400055)Tricyclics (400058)Zolpidem	Confirm (order): (4000126)Acetaminophen Confirmation (4000415)Amphetamine Class Confirmation (4000455)Barbiturate Class Confirmation (4000425)Benzodiazepine Class Confirmation (1000127)Buprenorphine and Opiate Analog Confirmation (4000465)Cannabinoid Class Confirmation (4000130)Carisoprodol Confirmation (4000435)Cocaine Class Confirmation (4000129)Fentanyl and Analog Confirmation (4000132)Gabapentin Confirmation (4000131)Methadone Confirmation (4000140)Methylphenidate Confirmation (4000445)Opiate Class Confirmation (4000146)Phencyclidine Confirmation (4000136)Pregabalin Confirmation (4000133)Propoxyphene Confirmation (4000138)Pseudoepheherine Confirmation (4000135)Tapentadol Confirmation (4000137)Tramadol Confirmation (4000142)Trazodone Confirmation (4000128)Tricyclic Confirmation (4000144)Zolpidem Confirmation	59384-8	4000070	Acetaminophen	Acetaminophen Screen	N/A	Alpha		
		N/A	4000090	Acetam Scr Int	Acetaminophen Screen Interp	N/A	Alpha		
		19343-3	4000051	Amphetamines	Amphetamines Screen	N/A	Alpha		
		N/A	4000091	Amphet Scr Int	Amphetamines Screen Interp	N/A	Alpha		
		19270-8	4000052	Barbiturates	Barbiturates Screen	N/A	Alpha		
		N/A	4000092	Barbit Scr Int	Barbiturates Screen Interp	N/A	Alpha		
		14316-4	4000053	Benzodiazepines	Benzodiazepines Screen	N/A	Alpha		
		N/A	4000093	Benzo Scr Int	Benzodiazepines Screen Interp	N/A	Alpha		
		58359-1	4000060	Buprenorphine	Buprenorphine and Opiate Analog Screen	N/A	Alpha		
		N/A	4000061	Bupren Scr Int	Buprenorphine and Opiate Analog Screen Inte	N/A	Alpha		
		18282-4	4000081	Cannabinoids	Cannabinoids Screen	N/A	Alpha		
		N/A	4000094	Cannabi Scr Int	Cannabinoids Screen Interp	N/A	Alpha		
		58370-8	4000056	Carisoprodol	Carisoprodol Screen	N/A	Alpha		
		N/A	4000095	Cariso Scr Int	Carisoprodol Screen Interp	N/A	Alpha		
		42241-0	4000080	Cocaine	Cocaine Screen	N/A	Alpha		
		N/A	4000096	Cocaine Scr Int	Cocaine Screen Interp	N/A	Alpha		
		58379-9	4000062	Fentanyl	Fentanyl and Analog Screen	N/A	Alpha		
		N/A	4000063	Fentanyl ScrInt	Fentanyl and Analog Screen Interp	N/A	Alpha		
		68439-9	4000059	Gabapentin	Gabapentin Screen	N/A	Alpha		
		N/A	4000097	Gabapen Scr Int	Gabapentin Screen Interp	N/A	Alpha		
		19550-3	4000064	Methadone	Methadone Screen	N/A	Alpha		
		N/A	4000065	Methadone Scr Int	Methadone Screen Interp	N/A	Alpha		
		19577-6	4000084	Methylphenidate	Methylphenidate Screen	N/A	Alpha		
		N/A	4000098	Methylp Scr Int	Methylphenidate Screen Interp	N/A	Alpha		
		19295-5	4000054	Opiates	Opiates Screen	N/A	Alpha		
		N/A	4000099	Opiate Scr Int	Opiate Screen Interp	N/A	Alpha		
		19659-2	4000079	Phencyclidine	Phencyclidine Screen	N/A	Alpha		
		N/A	4000100	Phencyc Scr Int	Phencyclidine Screen Interp	N/A	Alpha		
		64137-3	4000057	Pregabalin	Pregabalin Screen	N/A	Alpha		
		N/A	4000101	Pregab Scr Int	Pregabalin Screen Interp	N/A	Alpha		
		19429-0	4000066	Propoxyphene	Propoxyphene Screen	N/A	Alpha		
		N/A	4000067	Propox Scr Int	Propoxyphene Screen Interp	N/A	Alpha		
		4003-0	4000083	Pseudoephedrine	Pseudoephedrine Screen	N/A	Alpha		
		N/A	4000102	Pseudo Scr Int	Pseudoephedrine Screen Interp	N/A	Alpha		
		72485-6	4000068	Tapentadol	Tapentadol Screen	N/A	Alpha		
		N/A	4000069	Tapenta Scr Int	Tapentadol Screen Interp	N/A	Alpha		
		19710-3	4000048	Tramadol	Tramadol Screen	N/A	Alpha		
		N/A	4000049	Tramadol ScrInt	Tramadol Screen Interp	N/A	Alpha		
		4065-9	4000082	Trazodone	Trazodone Screen	N/A	Alpha		
		N/A	4000103	Trazad Scr Int	Trazodone Screen Interp	N/A	Alpha		
19312-8	4000055	Tricyclics	Tricyclics Screen	N/A	Alpha				
N/A	4000104	Tricyc Scr Int	Tricyclics Screen Interp	N/A	Alpha				
33338-5	4000058	Zolpidem	Zolpidem Screen	N/A	Alpha				
N/A	4000105	Zolpid Scr Int	Zolpidem Screen Interp	N/A	Alpha				
54247-2	4000075	Drug Screen Ur	Pain Panel Interp	N/A	Alpha				

New Test Effective : 11/27/2017										
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
4000315	.Pain Panel General Urine	.PAIN PNLG	Group							
CPT: 80307 HCPCS: G0479 Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months A "Detected" result on the following screen components will reflex the indicated confirmation order: Screen (assay): (4000070)Acetaminophen (4000051)Amphetamines (4000052)Barbiturates (4000053)Benzodiazepines (4000060)Buprenorphine (4000056)Carisoprodol (4000062)Fentanyl (4000059)Gabapentin (4000064)Methadone (4000054)Opiates (4000057)Pregabalin (4000066)Propoxyphene (4000068)Tapentadol (4000048)Tramadol (4000055)Tricyclics (4000058)Zolpidem	Confirmation (order):									
			(4000126)Acetaminophen Confirmation		59384-8	4000070	Acetaminophen	Acetaminophen Screen	N/A	Alpha
			(4000415)Amphetamine Class Confirmation		N/A	4000090	Acetam Scr Int	Acetaminophen Screen Interp	N/A	Alpha
			(4000455)Barbiturate Class Confirmation		19343-3	4000051	Amphetamines	Amphetamines Screen	N/A	Alpha
			(4000425)Benzodiazepine Class Confirmation		N/A	4000091	Amphet Scr Int	Amphetamines Screen Interp	N/A	Alpha
			(1000127)Buprenorphine and Opiate Analog Confirmation		19270-8	4000052	Barbiturates	Barbiturates Screen	N/A	Alpha
			(4000130)Carisoprodol Confirmation		N/A	4000092	Barbit Scr Int	Barbiturates Screen Interp	N/A	Alpha
			(4000129)Fentanyl and Analog Confirmation		14316-4	4000053	Benzodiazepines	Benzodiazepines Screen	N/A	Alpha
			(4000132)Gabapentin Confirmation		N/A	4000093	Benzo Scr Int	Benzodiazepines Screen Interp	N/A	Alpha
			(4000131)Methadone Confirmation		58359-1	4000060	Buprenorphine	Buprenorphine and Opiate Analog Screen	N/A	Alpha
			(4000445)Opiate Class Confirmation		N/A	4000061	Bupren Scr Int	Buprenorphine and Opiate Analog Screen Inte	N/A	Alpha
			(4000136)Pregabalin Confirmation		58370-8	4000056	Carisoprodol	Carisoprodol Screen	N/A	Alpha
			(4000133)Propoxyphene Confirmation		N/A	4000095	Cariso Scr Int	Carisoprodol Screen Interp	N/A	Alpha
			(4000135)Tapentadol Confirmation		58379-9	4000062	Fentanyl	Fentanyl and Analog Screen	N/A	Alpha
			(4000128)Tricyclic Confirmation		N/A	4000063	Fentanyl ScrInt	Fentanyl and Analog Screen Interp	N/A	Alpha
			(4000144)Zolpidem Confirmation		68439-9	4000059	Gabapentin	Gabapentin Screen	N/A	Alpha
					N/A	4000097	Gabapen Scr Int	Gabapentin Screen Interp	N/A	Alpha
					19550-3	4000064	Methadone	Methadone Screen	N/A	Alpha
					N/A	4000065	Methadone Scr Int	Methadone Screen Interp	N/A	Alpha
					19295-5	4000054	Opiates	Opiates Screen	N/A	Alpha
					N/A	4000099	Opiate Scr Int	Opiate Screen Interp	N/A	Alpha
					64137-3	4000057	Pregabalin	Pregabalin Screen	N/A	Alpha
					N/A	4000101	Pregab Scr Int	Pregabalin Screen Interp	N/A	Alpha
					19429-0	4000066	Propoxyphene	Propoxyphene Screen	N/A	Alpha
					N/A	4000067	Propox Scr Int	Propoxyphene Screen Interp	N/A	Alpha
					72485-6	4000068	Tapentadol	Tapentadol Screen	N/A	Alpha
					N/A	4000069	Tapenta Scr Int	Tapentadol Screen Interp	N/A	Alpha
					19710-3	4000048	Tramadol	Tramadol Screen	N/A	Alpha
					N/A	4000049	Tramadol ScrInt	Tramadol Screen Interp	N/A	Alpha
					19312-8	4000055	Tricyclics	Tricyclics Screen	N/A	Alpha
				N/A	4000104	Tricyc Scr Int	Tricyclics Screen Interp	N/A	Alpha	
				33338-5	4000058	Zolpidem	Zolpidem Screen	N/A	Alpha	
				N/A	4000105	Zolpid Scr Int	Zolpidem Screen Interp	N/A	Alpha	
				54247-2	4000075	Drug Screen Ur	Pain Panel Interp	N/A	Alpha	
New Test Effective : 11/27/2017										
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
4000126	Acetaminophen Confirmation Urine	ACETAM CNF	Group							
HCPCS: G0480 (1 Class) Class Information Below .Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.	Confirmation (order):									
					13622-6	4000246	Acetaminophen	Acetaminophen Confirmation	<250 ng/mL	XXXXX
					N/A	4000247	Acetaminoph Int	Acetaminophen Confirmation Interp	N/A	Alpha

New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000415	Amphetamine Class Confirmation Urine	AMPH CONF	Group						
HCPCS: G0480 (2 Classes) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			16234-7	4000152	Amphetamine	Amphetamine Confirmation	<50 ng/mL	XXXXX
				N/A	4000153	Amphet Int	Amphetamine Confirmation Interp	N/A	Alpha
				16235-4	4000154	Methamphetamin	Methamphetamine Confirmation	<100 ng/mL	XXXXXX
				N/A	4000155	Methamphet Int	Methamphetamine Confirmation Interp	N/A	Alpha
				18358-2	4000258	MDMA	MDMA Confirmation	<50 ng/mL	XXXXX
				N/A	4000259	MDMA Int	MDMA Confirmation Interp	N/A	Alpha
				20557-5	4000266	Phentermine	Phentermine Confirmation	<20 ng/mL	XXXXX
				N/A	4000267	Phentermine Int	Phentermine Confirmation Interp	N/A	Alpha
.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.									
New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000455	Barbiturate Class Confirmation Urine	BARBITCONF	Group						
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			16236-2	4000156	Butabarbitol	Butabarbitol Confirmation	<200 ng/mL	XXXXXX
				N/A	4000157	Butbarbitol Int	Butabarbitol Confirmation Interp	N/A	Alpha
				16237-0	4000158	Butalbital	Butalbital Confirmation	<200 ng/mL	XXXXXX
				N/A	4000159	Butalbital Int	Butalbital Confirmation Interp	N/A	Alpha
				16240-4	4000160	Pentobarbital	Amo/Pentobarbital Confirmation	<200 ng/mL	XXXXXX
				N/A	4000161	Pentobarb Int	Amo/Pentobarbital Confirmation Interp	N/A	Alpha
				16241-2	4000162	Phenobarbital	Phenobarbital Confirmation	<200 ng/mL	XXXXXX
				N/A	4000163	Phenobarbit Int	Phenobarbital Confirmation Interp	N/A	Alpha
				16238-8	4000290	Secobarbital	Secobarbital Confirmation	<200 ng/mL	XXXXXX
				N/A	4000291	Secobarb Int	Secobarbital Confirmation Interp	N/A	Alpha
.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.									

New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000425	Benzodiazepine Class Confirmation Urine	BENZO CONF	Group						
HCPCS: G0480 (1 Class) Class Information Below .Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.	Collect: Random Urine			16348-5	4000164	A-OH-alprazolam	A-OH-alprazolam Confirmation	<50 ng/mL	XXXXX
	Instructions: Sterile Container or Monovette 10mL (Min: 1mL)			N/A	4000165	A-OH-alpraz Int	A-OH-alprazolam Confirmation Interp	N/A	Alpha
	Storage/Transport Temp: Refrigerated			49876-6	4000166	A-OH-triazolam	A-OH-triazolam Confirmation	<50 ng/mL	XXXXX
	Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles			N/A	4000167	A-OH-triaz Int	A-OH-triazolam Confirmation Interp	N/A	Alpha
	Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			51777-1	4000168	7A-flunitrazep	7A-flunitrazepam Confirmation	<50 ng/mL	XXXXX
				N/A	4000169	7A-flunitra Int	7A-flunitrazepam Confirmation Interp	N/A	Alpha
				51776-3	4000170	7A-clonazepam	7A-clonazepam Confirmation	<50 ng/mL	XXXXX
				N/A	4000171	7A-clonaz Int	7A-clonazepam Confirmation Interp	N/A	Alpha
				59615-5	4000172	Alprazolam	Alprazolam Confirmation	<50 ng/mL	XXXXX
				N/A	4000173	Alprazolam Int	Alprazolam Confirmation Interp	N/A	Alpha
				16229-7	4000174	Clonazepam	Clonazepam Confirmation	<50 ng/mL	XXXXX
				N/A	4000175	Clonazepam Int	Clonazepam Confirmation Interp	N/A	Alpha
				16230-5	4000176	Desalkylfluraz	Desalkylflurazepam Confirmation	<100 ng/mL	XXXXXX
				N/A	4000177	Desalkylflura Int	Desalkylflurazepam Confirmation Interp	N/A	Alpha
				16227-1	4000178	Diazepam	Diazepam Confirmation	<100 ng/mL	XXXXXX
				N/A	4000179	Diazepam Int	Diazepam Confirmation Interp	N/A	Alpha
				20528-6	4000180	Flunitrazepam	Flunitrazepam Confirmation	<50 ng/mL	XXXXX
				N/A	4000181	Flunitrazep Int	Flunitrazepam Confirmation Interp	N/A	Alpha
				16231-3	4000182	Flurazepam	Flurazepam Confirmation	<50 ng/mL	XXXXX
				N/A	4000183	Flurazepam Int	Flurazepam Confirmation Interp	N/A	Alpha
				17088-6	4000184	Lorazepam	Lorazepam Confirmation	<50 ng/mL	XXXXX
				N/A	4000185	Lorazepam Int	Lorazepam Confirmation Interp	N/A	Alpha
				16228-9	4000186	Nordiazepam	Nordiazepam Confirmation	<25 ng/mL	XXXXX
				N/A	4000187	Nordiazepam Int	Nordiazepam Confirmation Interp	N/A	Alpha
				16201-6	4000188	Oxazepam	Oxazepam Confirmation	<50 ng/mL	XXXXX
				N/A	4000189	Oxazepam Int	Oxazepam Confirmation Interp	N/A	Alpha
				20559-1	4000190	Temazepam	Temazepam Confirmation	<50 ng/mL	XXXXX
				N/A	4000191	Temazepam Int	Temazepam Confirmation Interp	N/A	Alpha
			16232-1	4000192	Triazolam	Triazolam Confirmation	<25 ng/mL	XXXXX	
			N/A	4000193	Triazolam Int	Triazolam Confirmation Interp	N/A	Alpha	

New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000127	Buprenorphine and Opiate Analog Confirmation Urine	Bup Ana Conf	Group						
HCPCS: G0480 (2 Classes) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			49752-9	4000210	Buprenorphine	Buprenorphine Confirmation	<10 ng/mL	XXXXX
				N/A	4000211	Buprenor Int	Buprenorphine Confirmation Interp	N/A	Alpha
				49751-1	4000228	Norbuprenorphin	Norbuprenorphine Confirmation	<10 ng/mL	XXXXX
				N/A	4000229	Norbupren Int	Norbuprenorphine Confirmation Interp	N/A	Alpha
				16253-7	4000220	Meperidine	Meperidine Confirmation	<20 ng/mL	XXXXX
				N/A	4000221	Meperdin Int	Meperidine Confirmation Interp	N/A	Alpha
				58389-8	4000286	Normeperidine	Normeperidine Confirmation	<25 ng/mL	XXXXX
				N/A	4000287	Normeperdin Int	Normeperidine Confirmation Interp	N/A	Alpha
				42618-9	4000224	Naloxone	Naloxone Confirmation	<50 ng/mL	XXXXX
				N/A	4000225	Naloxone Int	Naloxone Confirmation Interp	N/A	Alpha
19610-5	4000226	Naltrexone	Naltrexone Confirmation	<50 ng/mL	XXXXX				
N/A	4000227	Naltrexone Int	Naltrexone Confirmation Interp	N/A	Alpha				
.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.									
New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000465	Cannabinoid Class Confirmation Urine	CANNABCONF	Group						
HCPCS: G0480 (2 Classes) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			20521-1	4000194	THC-11-Nor-9	THC-11-Nor-9 Confirmation	<25 ng/mL	XXXXX
				N/A	4000195	THC-11-Nor-9 Int	THC-11-Nor-9 Confirmation Interp	N/A	Alpha
				72461-7	4000196	JWH-018	JWH-018 Confirmation	<10 ng/mL	XXXXX
				N/A	4000197	JWH-018 Int	JWH-018 Confirmation Interp	N/A	Alpha
				72778-4	4000198	JWH-073	JWH-073 Confirmation	<10 ng/mL	XXXXX
				N/A	4000199	JWH-073 Int	JWH-073 Confirmation Interp	N/A	Alpha
.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.									
New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000130	Carisoprodol Confirmation Urine	CARISO CNF	Group						
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			58427-6	4000252	Carisoprodol	Carisoprodol Confirmation	<50 ng/mL	XXXXX
				N/A	4000253	Carisopro Int	Carisoprodol Confirmation Interp	N/A	Alpha
				58374-0	4000260	Meprobamate	Meprobamate Confirmation	<50 ng/mL	XXXXX
				N/A	4000261	Meprobamate Int	Meprobamate Confirmation Interp	N/A	Alpha
.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.									

New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000435	Cocaine Class Confirmation Urine	COCAI CONF	Group						
Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Class Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Information Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months Below				16250-3	4000202	Cocain Metab	Cocaine Confirmation	<50 ng/mL	XXXXX
				N/A	4000203	Cocai Metab Int	Cocaine Confirmation Interp	N/A	Alpha
.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.									
New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000129	Fentanyl and Analog Confirmation Urine	Fentanyl Conf	Group						
Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Class Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Information Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months Below				58381-5	4000214	Fentanyl	Fentanyl Confirmation	<2 ng/mL	XXXX
				N/A	4000215	Fentanyl Int	Fentanyl Confirmation Interp	N/A	Alpha
				58383-1	4000230	Norfentanyl	Norfentanyl Confirmation	<5 ng/mL	XXXX
				N/A	4000231	Norfentanyl Int	Norfentanyl Confirmation Interp	N/A	Alpha
				9396-3	4000274	Sufentanil	Sufentanil Confirmation	<5 ng/mL	XXXX
				N/A	4000275	Sufentanil Int	Sufentanil Confirmation Interp	N/A	Alpha
.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.									
New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000132	Gabapentin Confirmation Urine	GABAPE CNF	Group						
Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Class Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Information Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months Below				72810-5	4000254	Gabapentin	Gabapentin Confirmation	<500 ng/mL	XXXXXX
				N/A	4000255	Gabapentin Int	Gabapentin Confirmation Interp	N/A	Alpha
.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.									

New Test Effective : 11/27/2017											
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type		
4000140	Methylphenidate Confirmation Urine	MPHENI CNF	Group								
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			20548-4	4000262	Methylphenidate	Methylphenidate Confirmation	<50 ng/mL	XXXXX		
				N/A	2000263	Methylphen Int	Methylphenidate Confirmation Interp	N/A	Alpha		
				72790-9	4000272	Ritalinic Acid	Ritalinic Acid Confirmation	<50 ng/mL	XXXXX		
				N/A	4000273	Rital Acid Int	Ritalinic Acid Confirmation Interp	N/A	Alpha		
.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.											
New Test Effective : 11/27/2017											
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type		
4000131	Methadone Confirmation Urine	Methadone Conf	Group								
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			16246-1	4000206	Methadone	Methadone Confirmation	<50 ng/mL	XXXXX		
				N/A	4000207	Methadone Int	Methadone Confirmation Interp	N/A	Alpha		
				58429-2	4000204	EDDP	EDDP Confirmation	<100 ng/mL	XXXXXX		
				N/A	4000205	EDDP Int	EDDP Confirmation Interp	N/A	Alpha		
.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.											
New Test Effective : 11/27/2017											
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type		
4000108	Opiate Class Confirmation Urine	OPIATE CONF	Group								
HCPCS: G0480 (3 Classes) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			19593-3	4000208	Heroin Metab	Heroin Confirmation	<10 ng/mL	XXXXX		
				N/A	4000209	Heroin Met Int	Heroin Confirmation Interp	N/A	Alpha		
				16250-3	4000212	Codeine	Codeine Confirmation	<50 ng/mL	XXXXX		
				N/A	4000213	Codeine Int	Codeine Confirmation Interp	N/A	Alpha		
				16252-9	4000216	Hydrocodone	Hydrocodone Confirmation	<50 ng/mL	XXXXX		
				N/A	4000217	Hydrocodone Int	Hydrocodone Confirmation Interp	N/A	Alpha		
				16998-7	4000218	Hydromorphone	Hydromorphone Confirmation	<50 ng/mL	XXXXX		
				N/A	4000219	Hydromorp Int	Hydromorphone Confirmation Interp	N/A	Alpha		
				20550-0	4000222	Morphine	Morphine Confirmation	<50 ng/mL	XXXXX		
				N/A	4000223	Morphine Int	Morphine Confirmation Interp	N/A	Alpha		
				61421-4	4000232	Norhydrocodone	Norhydrocodone Confirmation	<50 ng/mL	XXXXX		
				N/A	4000233	Norhydrocod Int	Norhydrocodone Confirmation Interp	N/A	Alpha		
				61426-5	4000236	Noroxycodone	Noroxycodone Confirmation	<50 ng/mL	XXXXX		
				N/A	4000237	Noroxycodon Int	Noroxycodone Confirmation Interp	N/A	Alpha		
				16249-5	4000238	Oxycodone	Oxycodone Confirmation	<50 ng/mL	XXXXX		
				N/A	4000239	Oxycodone Int	Oxycodone Confirmation Interp	N/A	Alpha		
				17395-5	4000240	Oxymorphone	Oxymorphone Confirmation	<50 ng/mL	XXXXX		
				N/A	4000241	Oxymorphone Int	Oxymorphone Confirmation Interp	N/A	Alpha		
		.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.									

New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000146	Phencyclidine Confirmation Urine	PCP CNF	Group	16254-5	4000288	Phencyclidine	Phencyclidine Confirmation	<10 ng/mL	XXXXX
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			N/A	4000289	Phencyc Int	Phencyclidine Confirmation Interp	N/A	Alpha
				.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.					
New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000136	Pregabalin Confirmation Urine	PREGAB CNF	Group	64125-8	4000268	Pregabalin	Pregabalin Confirmation	<50 ng/mL	XXXXXX
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			N/A	4000269	Pregabalin Int	Pregabalin Confirmation Interp	N/A	Alpha
				.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.					
New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000133	Propoxyphene Confirmation Urine	Propoxy Conf	Group	16242-0	4000244	Propoxyphene	Propoxyphene Confirmation	<50 ng/mL	XXXXX
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			N/A	4000245	Propoxyphen Int	Propoxyphene Confirmation Interp	N/A	Alpha
				17377-3	4000242	Norpropoxyphene	Norpropoxyphene Confirmation	<50 ng/mL	XXXXX
				N/A	4000243	Norpropoxy Int	Norpropoxyphene Confirmation Interp	N/A	Alpha
.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.									
New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000138	Pseudoephedrine Confirmation Urine	PSEUDO CNF	Group	4004-8	4000270	Pseudoephedrine	Pseudoephedrine Confirmation	<20 ng/mL	XXXXX
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			N/A	4000271	Pseudoephed Int	Pseudoephedrine Confirmation Interp	N/A	Alpha
				.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.					
Only qualifies as 1 billable "class" if ordered alone. If ordered in conjunction with the Amphetamine Class Confirmation Urine (4000415) the two orderables combine into 2 total billable classes not 3.									

New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000135	Tapentadol Confirmation Urine	Tapentadol Conf	Group						
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			65807-0	4000294	Tapentadol	Tapentadol Confirmation	<50 ng/mL	XXXXX
				N/A	4000295	Tapentadol Int	Tapentadol Confirmation Interp	N/A	Alpha
				65808-8	4000296	NDesMethTap	N-DesMethylTapentadol Confirmation	<50 ng/mL	XXXXX
				N/A	4000297	NDesMethTap Int	N-DesMethylTapentadol Confirmation Interp	N/A	Alpha
.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.									
New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000137	Tramadol Confirmation Urine	Tramadol Conf	Group						
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			20561-7	4000234	Tramadol	Tramadol Confirmation	<50 ng/mL	XXXXX
				N/A	4000235	Tramadol Int	Tramadol Confirmation Interp	N/A	Alpha
				.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.					
New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000142	Trazodone Confirmation Urine	TRAZOD CNF	Group						
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine Instructions: Sterile Container or Monovette 10mL (Min: 1mL) Storage/Transport Temp: Refrigerated Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			59975-3	4000280	Trazodone	Trazodone Confirmation	<25 ng/mL	XXXXX
				N/A	4000281	Trazodone Int	Trazodone Confirmation Interp	N/A	Alpha
				.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.					

New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000128	Tricyclic Confirmation Urine	TRICYC CNF	Group	20515-3	4000248	Amitriptyline	Amitriptyline Confirmation	<50 ng/mL	XXXXX
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine			N/A	4000249	Amitript Int	Amitriptyline Confirmation Interp	N/A	Alpha
	Instructions: Sterile Container or Monovette 10mL (Min: 1mL)			61419-8	4000256	Imipramine	Imipramine Confirmation	<50 ng/mL	XXXXX
	Storage/Transport Temp: Refrigerated			N/A	4000257	Imipramine Int	Imipramine Confirmation Interp	N/A	Alpha
	Unacceptable Conditions: Tubes with additives, repeated freeze-thaw cycles			61428-9	4000264	Nortriptyline	Nortriptyline Confirmation	<50 ng/mL	XXXXX
	Stability: Ambient: 7 days, Refrigerated: 30 days, Frozen: 6 months			N/A	4000265	Nortript Int	Nortriptyline Confirmation Interp	N/A	Alpha

.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.

New Test Effective : 11/27/2017									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4000144	Zolpidem Confirmation Urine	ZOLPID CNF	Group	72770-1	4000282	Zolpidem	Zolpidem Confirmation	<20 ng/mL	XXXXX
HCPCS: G0480 (1 Class) Class Information Below	Collect: Random Urine			N/A	4000283	Zolpidem Int	Zolpidem Confirmation Interp	N/A	Alpha
	Instructions: Sterile Container or Monovette 10mL (Min: 1mL)								

.Prescriptions (4000106) and Creatinine Urine (3000750) orders will be reflexed if ordered as a stand-alone confirmation. Confirmation orders reflexing from Pain Panel General (4000314) or Pain Panel Comprehensive (4000334) orders will not generate additional .Prescription (4000106) or Creatinine Urine (3000750) orders.

Definitive drug class testing is billed as a combination of the total number of classes tested. The combinations are as follows:

Total number of classes tested:	Billing Code:
1-7 Total Classes Tested	G0480
8-14 Total Classes Tested	G0481
15-21 Total Classes Tested	G0482
22+ Total Classes Tested	G0483

INACTIVATED TESTS - Please update your EMR catalog with those appropriate to your practice

INACTIVE Effective : Immediately						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
4000134	Meprobamate Confirmation	MEPROB CNF	4000130	CARISO CNF	Carisoprodol Confirmation	

INACTIVE Effective : 11/27/2017						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
4000445	Opiate Class Confirmation	OPIATE CONF	4000127	Bup Ana Conf	Buprenorphine and Opiate Analog Confirmation Urine	
			4000129	Fentanyl Conf	Fentanyl and Analog Confirmation Urine	
			4000131	Methadone Conf	Methadone Confirmation Urine	
			4000108	OPIATE CONF	Opiate Class Confirmation Urine	
			4000133	Propoxy Conf	Propoxyphene Confirmation Urine	
			4000135	Tapentadol Conf	Tapentadol Confirmation Urine	
			4000137	Tramadol Conf	Tramadol Confirmation Urine	