

NEW TESTS - Please update your EMR catalog with those appropriate to your practice

New Test Effective : 08/19/2019									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6907385	Blastomyces Antigen Quant by EIA, Urine	Blasto Ag Ur	Group	N/A	6907387	Blasto Ag Ur	Blastomyces Antigen Urine	By Report	Alpha
CPT: 87449	Specimen: 1mL (0.5mL) urine in a sterile container. Unacceptable: Specimens preserved in boric acid Stability: Ambient -2 weeks; Refrigerated -2 weeks; Frozen -Indefinitely Set-up: Tue, Thu, Sat TAT: 1-4 days after set-up Method: Quantitative Enzyme Immunoassay			N/A	6907389	BlastoAgU Intrap	Blastomyces dermatitidis EIA Interp	By Report	Alpha

New Test Effective : 09/02/2019									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
3870659	Neutrophil Oxidative Index	Neutro Ox Index	Group	n/a	5568850	Particle Stim	Particle Stimulation	>= 65 NOI	Numeric XXX.X
CPT: 86352	Collect Monday through Thursday only! Collect 4 mL (1 ml) Whole Blood in a EDTA (Lavender Top) DO NOT Centrifuge or Refrigerate Deliver to laboratory (flow cytometry) ASAP! Specimen stability: 24hrs Room Temperature.			85374-7	5568875	Soluble Stim	Soluble Stimulation	>= 329 NOI	Numeric XXX.X
				n/a	5568876	NOI Interp	NOI Interpretation	Normal	Alpha
This will be reported out in General Laboratory without a pathology interpretation									
Specimen must be received in the lab section before 2pm and within 8hrs of collection.									

MODIFICATION TESTS - Please update your EMR catalog with those appropriate to your practice

Modify Test Effective : 08/19/19 - Please update your EMR catalog with those appropriate to your practice -									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
5570915	Anti Mullerian Hormone	AMH	Detail	38476-8	5570915	AMH	Anti Mullerian Hormone	(see below) ng/mL	Numeric XX.XXX
CPT: 83520	Test will now be performed in house. Specimen: Serum 0.5 mL(0.2 mL) Clot Activator SST Instructions: Allow serum to clot then Separate from cells ASAP. Unacceptable Conditions: Hemolyzed or lipemic specimens. Stability: Ambient 1 day, Refrigerated 6 days, Frozen > 6 days Testing Schedule: Mon, Wed, Fri								
							Age Range	Male	Female
							0 - 18 Years	n/a	n/a
							18 - 26 Years	n/a	1.02 - 14.63
							26 - 31 Years	n/a	0.69 - 13.39
							31 - 36 Years	n/a	0.36 - 10.07
							36 - 41 Years	n/a	0.18 - 5.68
							41 - 46 Years	n/a	0.01 - 2.99
							46 - 150 Years	n/a	n/a

Modify Test Effective : 08/19/2019 - Please update your EMR catalog with those appropriate to your practice -									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6906641	Zika Virus IgM	Zika Virus IgM	Group	80824-6	6906643	Zika Virus IgM	Zika Virus IgM	Negative	Alpha
	This test will reflex Zika IgM Conf if result is reported as "Presumptive Zika". See below for test build information No other test changes			66174-4	6906645	Pregnant Y/N	Pregnant Y/N	Prompt	Alpha
				NA	6906647	Zika Expose Y/N	Relevant Zika Exposure Y/N	Prompt	Alpha
				80824-6	6906649	Symptoms Y/N	Symptoms consistent with Zika Y/N	Prompt	Alpha

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6907209	Reflex Zika IgM Ab Capture (MAC) Confirmation	Zika IgM Conf	Detail	N/A	6907209	Zika IgM Conf	Zika IgM Ab Capture (MAC) Confirmation	Negative	Alpha
	This is the reflex test for Zika Virus IgM (6906641)								

Modify Test Effective : 08/19/2019 - Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6906191	Myositis Extended Panel	Myositis Panel	Group						
CPT: 83516x12 CPT changing from 83516x13, 86235x6 to 83516x12, 86235x6 86235x6 Test schedule change from Mon, Tues, Thurs to Sun-Sat One result component has been discontinued No other test changes									
			Discontinue	31625-7	5513216	SSA 52 (Ro) IgG	SSA 52 (Ro) ENA Ab IgG	AU/ml	Numeric XX
				53019-6	5513215	SSA 60 (Ro) IgG	SSA 60 (Ro) Ab IgG	AU/ml	Numeric XX
				33983-8	5513217	Ribonucleic P Ab	Ribonucleic P Ab	AU/ml	Numeric XX
				33571-1	5513220	JO-1 Auto Ab	JO-1 Autoantibody	AU/ml	Numeric XX
				33771-7	5513140	PL-12 Auto Ab	PL-12 Autoantibody	Negative	Alpha
				33772-5	5513130	PL-7 Auto Ab	PL-7 Autoantibody	Negative	Alpha
				45152-6	5513190	OJ Auto Ab	OJ Autoantibody	Negative	Alpha
				45149-2	5513180	EJ Auto Ab	EJ Autoantibody	Negative	Alpha
				33921-8	5513210	SRP Auto Ab	Signal Recognition Particle (SRP) Autoantibody	Negative	Alpha
				18484-6	5513170	Ku Auto Ab	Ku Autoantibody	Negative	Alpha
				81732-0	3806050	Pm-1 Ab	PM-SCL Antibody	Negative	Alpha
				NA	5513213	U2 sn RNP Ab	U2 sn RNP Ab	Negative	Alpha
				49963-2	6906193	Fibrillarin IgG	Fibrillarin (U3 RNP) Ab, IgG	Negative	Alpha
				18485-3	5513160	Mi-2 Auto Ab	Mi-2 Autoantibody	Negative	Alpha
				NA	5513211	P155/140 Ab	P155/140 Ab	Negative	Alpha
				82448-2	6906195	TIF1-gamma Ab	TIF1-gamma Antibody	Negative	Alpha
				82440-9	6906197	SAE1 Antibody	SAE1 (SUMO activating enzyme) Antibody	Negative	Alpha
				82424-3	6906199	MDA5 Antibody	MDA5 (CADM-140) Antibody	Negative	Alpha
				82426-8	6906201	NXP-2 Ab	NXP-2 (Nuclear matrix protein-2) Ab	Negative	Alpha
				48767-8	6906203	Myositis Pn Int	Myositis Panel Interpretative Data	Negative	Alpha

Modify Test Effective : 08/19/2019 - Please update your EMR catalog with those appropriate to your practice -

Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6906831	Blastomyces dermatitidis Ag Quant, EIA	Blasto Ag Quant	Group						
Urine is no longer an acceptable specimen for this test See Blastomyces Antigen Quant by EIA, Urine in Green New Test section for replacement build information - No other test changes				6313-1	6906833	Blasto Ag	Blastomyces Ag	ng/ml	Numeric XXX.XX
				N/A	6906835	Blasto Ag Sourc	Blastomyces Ag Source	ng/ml	Numeric XXX.XX

INACTIVATED TESTS - Please update your EMR catalog with those appropriate to your practice

INACTIVE Effective : 08/19/2019

Inactive Test Code	Inactive Test Name	Inactive Mnemonic	Replacement Information						
6905095	Pancreatitis (PRSS1) Sequencing	PRSS1 FGS	<table border="1"> <thead> <tr> <th>Test Code</th> <th>Mnemonic</th> <th>Test Name</th> </tr> </thead> <tbody> <tr> <td>n/a</td> <td>n/a</td> <td>n/a</td> </tr> </tbody> </table>	Test Code	Mnemonic	Test Name	n/a	n/a	n/a
Test Code	Mnemonic	Test Name							
n/a	n/a	n/a							

This test is being discontinued due to low utilization

INACTIVE Effective : 08/19/2019

Inactive Test Code	Inactive Test Name	Inactive Mnemonic	Replacement Information						
6905161	Sulfhemoglobin, Whole Blood	Sulf Hgb WB	<table border="1"> <thead> <tr> <th>Test Code</th> <th>Mnemonic</th> <th>Test Name</th> </tr> </thead> <tbody> <tr> <td>n/a</td> <td>n/a</td> <td>n/a</td> </tr> </tbody> </table>	Test Code	Mnemonic	Test Name	n/a	n/a	n/a
Test Code	Mnemonic	Test Name							
n/a	n/a	n/a							

This test is being discontinued due to low utilization

INACTIVE Effective : 08/19/2019

Inactive Test Code	Inactive Test Name	Inactive Mnemonic	Replacement Information						
801040	Hemoglobin Fluid	SRS HGB	<table border="1"> <thead> <tr> <th>Test Code</th> <th>Mnemonic</th> <th>Test Name</th> </tr> </thead> <tbody> <tr> <td>n/a</td> <td>n/a</td> <td>n/a</td> </tr> </tbody> </table>	Test Code	Mnemonic	Test Name	n/a	n/a	n/a
Test Code	Mnemonic	Test Name							
n/a	n/a	n/a							

This test is being discontinued due to low utilization

INACTIVE Effective : 08/19/2019						
Inactive	Inactive	Inactive	<i>Replacement Information</i>			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
801050	Hematocrit Body Fluid	SRS HCT	n/a	n/a	n/a	
This test is being discontinued due to low utilization						
INACTIVE Effective : 09/02/2019						
Inactive	Inactive	Inactive	<i>Replacement Information</i>			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
5569200	Neutrophil Oxidative Index (NOI)	CHEMILUMIN	3870659	Neutro Ox Index	Neutrophil Oxidative Index	
This test with Pathology Interpretation is being discontinued and replaced with General Laboratory test at the right. For Build information for the Replacement test please see the Green New Test Build section of this LiveWire.						