

NEW TESTS - Please update your EMR catalog with those appropriate to your practice									
New Test	Effective : 05/31/2021								
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
6907649	17-Ketosteroids, Urine	17Ketosteroid U	Group						
CPT: 83586	- Collect: 24-hour urine. Refrigerate during collection. - Specimen Prep: Transfer two 4mL aliquots from a well-mixed 24-hour urine collection into 2 transport tubes containing 20 mg Sulfamic Acid (Min 3ml/aliquot). Adequate refrigeration is the most important aspect of specimen preservation. - Transport Temperature: Refrigerated - Remarks: record total volume and collection time interval on transport tube and test request form. - Stability: Ambient - 4 hours; Refrigerated - 2 weeks; Frozen - 1month - Set-up: Mon, Wed, Fri TAT: 1-5 days after set-up - Method: Quant Spectrophotometry			95543	3873099	CreatUr-per vol	Creatinine, Urine-per volume	mg/dL	By Report
				95847	3810375	CreatUr-per 24h	Creatinine, Urine - per 24 hr	mg/d	By Report
				21038-5	6907643	17Ketos-per vol	17-Ketosteroids - per volume	mg/L	By Report
				6766-0	6907645	17Ketos-per 24h	17-Ketosteroids - per 24 hr	mg/d	By Report
				19153-6	N/A	Total Volume	Total Volume	mL	By Report
				30211-7	N/A	Hrs Collected	Hrs Collected	hr	By Report
				Prompt code => 3809175					
				Prompt code => 3000651					
MODIFICATION TESTS - Please update your EMR catalog with those appropriate to your practice									
Modify Test Effective : Immediately - Please update your EMR catalog with those appropriate to your practice -									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
4005950	Vancomycin Level Trough	VANCO TROU	Detail	4092-3	4005950	Vanco Tr	Vancomycin Trough	10.0 - 20.0 ug/mL	Numeric XXXX.X
CPT: 80202	Reference Range Changing from >10.0 ug/mL to 10.0-20.0 ug/mL No other changes								
Modify Test Effective : Immediately - Please update your EMR catalog with those appropriate to your practice -									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
5000855	Fetal Hemoglobin, Quant	HGB F QUANT	Detail	4576-5	105275	Fetal HGB Ref	Fetal HGB (Hemoglobinopathy Eval)	no change	no change
CPT: 83020	Change in Name from "Fetal Hemoglobin, HPLC" to "Fetal Hemoglobin, Quant" CPT code changes from 83021 to 83020. No other changes								
Modify Test Effective : 05/17/2021 - Please update your EMR catalog with those appropriate to your practice -									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
3803525	RBC Magnesium	MAG RBC	Detail	2597-3	3803525	RBC Mag	RBC Magesium	3.6-7.5 mg/dL	Numeric XXXXX.X
CPT: 83735	Change the reference range and UOM for this test from 1.5-3.1 mmol/L to 3.6-7.5 mg/dL No other changes								
Modify Test Effective : 05/17/2021 - Please update your EMR catalog with those appropriate to your practice -									
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type
5000010	Hemochromatosis (HFE) 3 Mutations	HEMOCHROMA	Group						
Yellow ACD Tubes (A or B solution) is now an acceptable alternate specimen type Unacceptable Specimens: Frozen specimens in glass tubes. Updated Frozen stability: 1 month. No other changes				21695-2	5000020	C282Y (Genotype)	C282Y Hemochromatosis Mutation	no change	no change
				21696-0	5000015	H63D (Genotype)	H63D Hemochromatosis Mutation	no change	no change
				38380-2	5000025	S65C (Genotype)	S65C Hemochromatosis Mutation	no change	no change
				34519-9	5000030	HFE Interpretation	Hemochromatosis Mutation Interpretation	no change	no change

Modify Test Effective : 05/17/2021 - Please update your EMR catalog with those appropriate to your practice -										
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
9108985	HLA Typing For Narcolepsy	HLA NARCO	Detail	43291-4	9109000	Narco HLA Intrap	HLA Narcolepsy Interpretation	no change	no change	
Yellow ACD Tubes (A or B solution) is now an acceptable alternate specimen type Unacceptable Specimens: Frozen specimens in glass tubes. Updated Frozen stability: 1 month. No other changes										
Modify Test Effective : 05/17/2021 - Please update your EMR catalog with those appropriate to your practice -										
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
1517000	Plasminogen Act Inhibitor, Activity Genotype	PLAS ACT G	Detail	N/A	3615985	PAI-1 Interp	PAI-1 Interpretation	no change	no change	
Yellow ACD Tubes (A or B solution) is now an acceptable alternate specimen type Unacceptable Specimens: Frozen specimens in glass tubes. Updated Frozen stability: 1 month. No other changes										
Modify Test Effective : 05/17/2021 - Please update your EMR catalog with those appropriate to your practice -										
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
6900555	Propoxyphene and Metabolite S/P Quant	Propoxyphene S/P	Group	3543-6	6900557	Propoxyphene	Propoxyphene	no change	no change	
				3870-3	6900559	Norpropoxyphene	Norpropoxyphene	no change	no change	
Gray top tube is no longer an acceptable specimen type Specimen volume changes from 1mL(0.5mL) to 2mL(0.7mL). Updated stability: Ambient - 2 weeks; Refrigerated - 2 weeks; Frozen - 3 months Updated test set-up: varies; updated TAT: 8 - 11 days No other changes										
Modify Test Effective : 5/31/2021 - Please update your EMR catalog with those appropriate to your practice -										
Test Code	Test Name	Mnemonic	Category/Type	LOINC	Result Code	Mnemonic	Result Name	Reference Range UOM	Result Type	
3609525	West Nile Virus Antibodies IgG and IgM	WEST NILE	Group	38997-3	3610075	W Nile IgG	West Nile Virus Antibody IgG			
CPT: 86789; 86788 Change in the Result codes for the IgG and IgM for West Nile antibodies and change in specimen stability IgG changes from 3610075 to 6904769 IgM changes from 3610125 to 6904771 Stability After separation from cells: Ambient 48 hours; Refrigerated 2 weeks; Frozen 1 year (avoid repeated freeze/thaw cycles) No other changes				Remove	38997-3	3610075	W Nile IgG	West Nile Virus Antibody IgG		
				Add New	38997-3	6904769	West Nile IgG	West Nile Virus Ab, IgG, Serum	0.00 -1.29 IV	Numeric XXXX.XX
				Remove	38166-5	3610125	W Nile IgM	West Nile Virus Antibody IgM		
				Add New	38166-5	6904771	West Nile IgM	West Nile Virus Ab, IgM, Serum	0.00 -0.89 IV	Numeric XXXX.XX
INACTIVATED TESTS - Please update your EMR catalog with those appropriate to your practice										
INACTIVE Effective : Immediately										
Inactive	Inactive	Inactive	Replacement Information							
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name					
5537425	Blackjack Oak Tree IgE	BLACKJACK	2935600	TREE OAK R	Red Oak Tree IgE	Contact RML for test build information				
INACTIVE Effective : Immediately										
Inactive	Inactive	Inactive	Replacement Information							
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name					
9107025	PCA3 - Prostate Cancer Biomarker	PCA3 URINE	N/A	N/A	N/A	There is no available replacement for this test				

INACTIVE Effective : Immediately						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
3487702	LDL Subclasses	LDL Subcls	Option 1 5520001	NMR Particle Ct	NMR Particle Count Only	
			Option 2 5520007	LIPID/NMR Part	Lipid Group with NMR Particle Count	
<p>This test is discontinued at the performing reference laboratory There are two options for replacement tests, the NMR Particle Count, or the Lipid Group with NMR Particle Count</p>			<p>Contact RML for test build information</p>			
INACTIVE Effective : 05/17/2021						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6904741	Alpha-1-Microglobulin, Urine	Alpha-1-Microgl Ur	N/A	N/A	N/A	
			<p>There is no available replacement for this test</p>			
INACTIVE Effective : 05/17/2021						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6904871	Rett Syndrome (MECP2), Full Gene Sequencing	Rett Syndrome	N/A	N/A	N/A	
			<p>There is no available replacement for this test</p>			
INACTIVE Effective : 05/17/2021						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6900535	Norovirus Antigen Stool	Norovirus Ag	6006467	Norovirus PCR	Norovirus Groups 1,2 PCR	
			<p>Contact RML for test build information</p>			
INACTIVE Effective : 05/17/2021						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
3488150	Alpha Thalassemia, 7 Deletions	Alpha Thal	6904751	Alpha-Globin	Alpha-Globin (HBA1 & HBA2) Deletion/Duplication	
			<p>Contact RML for test build information</p>			
INACTIVE Effective : 05/31/2021						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
5578350	Honey IgE	Honey	N/A	N/A	N/A	
			<p>There is no available replacement for this test</p>			
INACTIVE Effective : 05/31/2021						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
6907125	Allspice IgE	Allspice	5529725	Clove	Clove IgE	
			<p>Contact RML for test build information</p>			
INACTIVE Effective : 05/31/2021						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
5559760	Orange Roughy IgE	ORANGE RUORANGE RU	N/A	N/A	N/A	
			<p>There is no available replacement for this test</p>			
INACTIVE Effective : 05/31/2021						
Inactive	Inactive	Inactive	Replacement Information			
Test Code	Test Name	Mnemonic	Test Code	Mnemonic	Test Name	
3656050	17-Ketosteroids, 24-Hour Urine	KETOSTER17	6907649	17Ketosteroid U	17-Ketosteroids, Urine	
			<p>See New Test section for test build information</p>			