

Effective January 21st, 2013

NEW TESTS - Please update your EMR catalog with those appropriate to your practice

NEW TEST		Available 01/21/2013							
Order Test Code	Test Name	Mnemonic	Category/Type	Result Codes	Result Name	Mnemonic	UOM	Reference Range	Result
3805327	METANEPHRINES FRACTIONATED P	METANEPH P	Group Test	3804227	METANEPHRINE, PLASMA	XMETANEPH	pg/mL	Below	Numeric
				3804277	NORMMETANEPHRINES, PLASMA	XNORMETAN			
				3804237	METANEPHRINE TOTAL PLASMA	XTOT METAN			
CPT Code: 83835									
		<i>Specimen Collection:</i> 2mL (1.5mL) Plasma From EDTA Lavender - Spin Separate and Freeze ASAP! Specimen Stability: Room temperature: 4 hours, Refrigerated: 14 days, Frozen: 14 days							
		Metanephrine: < OR = 57 pg/mL Normetanephrine: < OR = 148 pg/mL Total Metanephrine: < OR = 205 pg/mL							

INACTIVE TESTS - Please update your EMR catalog with those appropriate to your practice

INACTIVE		Effective 01/21/2012					
Order Test Code	Test Name	Mnemonic	Replacement Code	Order Test Code	Test Name	Mnemonic	Reason
3804325	METANEPHRINES FRACTIONATED P	METANEPH P	Group Test	3805327	METANEPHRINES FRACTIONATED P	METANEPH P	Test Restructuring